

```
.. module:: bvp
   :platform: Unix, Windows
   :synopsis: This module provides various functions to handle BVP signals.

.. moduleauthor:: Filipe Canento
```

Modules

bspyclone.filt	os	bspyclone.plux	bspyclone.tools
glob	bspyclone.peakd	scipy	unittest
numpy	pylab	bspyclone.sync	

Functions

features(Signal=None, SamplingRate=1000.0, Filter={})
 Retrieves relevant BVP signal features.

Kwargs:

- Signal (array): input signal.
- SamplingRate (float): sampling frequency (Hz).
- Filter (dict): filter parameters.

Kwrvals:

- Signal (array): output filtered signal (see notes 1).
- Amplitude (array): signal pulses amplitudes (in the units of the input signal). (TODO)
- Onset (array): indexes (or instants in seconds, see notes 2.b) of the pulses onsets.
- Peak (array): indexes (or instants in seconds, see notes 2.b) of the pulses peaks. (TODO)
- DicroticNotch (array): indexes (or instants in seconds, see notes 2.b) of the pulses dicrotic notchs. (TODO)
- IBI (array): Inter-Beat Intervals in msec (see notes 2.a).
- HR (array): Instantaneous Heart Rates in b.p.m. (see notes 2.a).
- mean (float): mean
- std (float): standard deviation
- var (float): variance
- skew (ndarray): skewness
- kurtosis (array): kurtosis
- ad (float): absolute deviation

Configurable fields:
 {"name": "bvp.features", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter"], "outputs": ["Signal", "Amplitude", "Onset", "Peak", "DicroticNotch", "IBI", "HR", "mean", "std", "var", "skew", "kurtosis", "ad"]}

See Also:

- [filt](#)
- [pulse](#)
- [tls.statsf](#)

Notes:

- 1 - If a filter is given as a parameter, then the returned keyworded values dict has a 'Signal' key.
- 2 - If the sampling rate is defined, then:
 - a) the returned keyworded values dict has keys 'IBI' and 'HR'.
 - b) keys 'onset', 'peak', and 'DicroticNotch' are converted to instants of occurrence in seconds.

Example:

```
bvp = ...
SamplingRate = ...
res = pulse(Signal=bvp, SamplingRate=SamplingRate)
```

References:

```
.. [1]
```

filt(Signal=None, SamplingRate=1000.0, UpperCutoff=8.0, LowerCutoff=1.0, Order=4.0)
 Filters an input BVP signal.

If only input signal is provide, it returns the filtered signal assuming a 1000Hz sampling frequency and the default filter parameters: low-pass filter with cutoff frequency of 8Hz followed by a high-pass filter with cutoff frequency of 1Hz.

Kwargs:

- Signal (array): input signal.

SamplingRate (float): sampling frequency (Hz).
UpperCutoff (float): Low-pass filter cutoff frequency (Hz).
LowerCutoff (float): High-pass filter cutoff frequency (Hz).
Order (int): Filter order.

Kwrvals:
Signal (array): output filtered signal.

Configurable fields:{"name": "bvp.filt", "config": {"UpperCutoff": "8.", "SamplingRate": "1000.", "LowerCutoff": "1.", "Order": "4."}, "input"

See Also:
flt.zpdfnr

Notes:

Example:

References:
.. [1]

onset(Signal=None, SamplingRate=1000.0)

Determines the onsets of the BVP signal pulses.
Skips very corrupted signal parts.

Kwargs:
Signal (array): input signal.
SamplingRate (float): sampling frequency (Hz).

Kwrvals:
Onset (array):

Configurable fields:{"name": "bvp.onset", "config": {"SamplingRate": "1000."}, "inputs": ["Signal"], "outputs": ["Onset"]}

See Also:

Notes:

Example:

References:
.. [1]

pulse(Signal=None, SamplingRate=1000.0, Filter={})

Determines BVP signal pulse information.

Kwargs:
Signal (array): input signal.
SamplingRate (float): sampling frequency (Hz).
Filter (dict): filter parameters.

Kwrvals:
Signal (array): output filtered signal (see notes 1).
Amplitude (array): signal pulses amplitudes (in the units of the input signal). (TODO)
Onset (array): indexes (or instants in seconds, see notes 2.b) of the pulses onsets.
Peak (array): indexes (or instants in seconds, see notes 2.b) of the pulses peaks. (TODO)
DicroticNotch (array): indexes (or instants in seconds, see notes 2.b) of the pulses dicrotic notchs. (TODO)
IBI (array): Inter-Beat Intervals in msec (see notes 2.a).
HR (array): Instantaneous Heart Rates in b.p.m. (see notes 2.a).

Configurable fields:{"name": "bvp.pulse", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter"], "outputs": ["Signal", "Amplitu"

See Also:
filt
onset

Notes:
1 - If a filter is given as a parameter, then the returned keyworded values dict has a 'Signal' key.
2 - If the sampling rate is defined, then:
a) the returned keyworded values dict has keys 'IBI' and 'HR'.
b) TODO: keys 'onset', 'peak', and 'DicroticNotch' are converted to instants of occurrence in seconds.

Example:
bvp = ...
SamplingRate = ...
res = [pulse](#)(Signal=bvp, SamplingRate=SamplingRate)

References:
.. [1]



```
.. module:: biomesh
    :platform: Unix, Windows
    :synopsis: This module provides an API to use the BioMESH specification. See the specification at http://camoes.lx.it.pt/MediaWiki/index.php/Database\_S

.. moduleauthor:: Carlos Carreiras
```

Modules

bspyclone.database.h5db	pymongo	string
numpy	re	sympy
os	bspyclone.database.syncdb	warnings

Classes

[biomesh](#)
[dataContainer](#)
[dataSelector](#)
[experiments](#)
[records](#)
[subjects](#)

class biomesh

Class to operate on a BioMESH DB database.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
 .. [1]

Methods defined here:

__init__(self, dbName=None, host='localhost', port=27017, dstPath='~', srvPath='/BioMESH', sync=False, altSync=False)
 Establish a connection to a BioMESH DB server. If the database (DB) does not exist, one is created with the necessary basic structures.

Kwargs:

dbName (str): Name of the DB to connect to.

host (str): Network address of the MongoDB server. Default: 'localhost'.

port (int): Port the MongoDB server is listening on. Default: 27017.

dstPath (str): Path to store the HDF5 files. Default: '~'.

srvPath (str): Path to store the HDF5 files on the remote server. Default: '/BioMESH'.

sync (bool): Flag to perform synchronization with remote server. Default: True.

altSync (bool): New experimental synchronization framework. Default: False.

Kwrvals:

See Also:

Notes:

If the 'sync' flag is True, the connection to the database may take longer to establish due to the synchronization step. This is es

Example:

db = [biomesh](#)(dbName='biomesh_tst', host='193.136.222.234', port=27017, dstPath='~/tmp/[biomesh](#)', srvPath='/biomesh_tst', sync=False)

References:
 .. [1]

close(self)

Close the connection to the database.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:
`db.close()`

References:
.. [1]

drop(self, flag=True)

Function to remove the database from the server. The HDF5 files are preserved.

Kwargs:
flag (bool): Flag to override user confirmation about removal. Set to False to avoid user prompt. Default: True.

Kwrvals:

See Also:

Notes:

Example:
`db.drop()`

References:
.. [1]

expsInSub(self, subjectId)

List all the [experiments](#) in the database belonging to a given subject.

Kwargs:
subjectId (int): The ID of the subject.

Kwrvals:
nameList (list): List with the names of the [experiments](#) in the database.

See Also:

Notes:

Example:
`explist = db.expsInSub(0)['nameList']`

References:
.. [1]

h5Add(self, filePath=None)

Procedure to add already created HDF5 files ([records](#)) to MongoDB.

Kwargs:
filePath (str): Location of the file to add.

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

metaH5Add(self, filePath=None)

Procedure to add an already created meta HDF5 file ([experiments](#) and [subjects](#)) to MongoDB.

Kwargs:
filePath (str): Location of the file to add.

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

subsInExp(self, experimentName)

List all the [subjects](#) in the database belonging to a given experiment.

Kwargs:
 experimentName (str): The name of the experiment.

Kwrvals:
 idList (list): List with the IDs of the [subjects](#).

See Also:

Notes:

Example:
 subList = db.[subsInExp](#)('unknown')['idList']

References:
 .. [1]

syncNow(self)

Function to perform the synchronization with the remote server.

Kwargs:
 None

Kwrvals:
 None

See Also:

Notes:

Example:
 db.sync()

References:
 .. [1]

class dataContainer

Data container for signals and events.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
 .. [1]

Methods defined here:

__init__(self, case=None, mdata={}, signal=[], timeStamps=[], values=[])
 Initialize the container according to case (signal or events).

Kwargs:
 case (str): Case ('signals' or 'events') to store.

Kwrvals:

See Also:

Notes:

Example:

References:
 .. [1]

__str__(self)
 str operator.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

class **dataSelector**

Wrapper for operator overloading to [records](#).

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

Methods defined here:

__getitem__(self, key)
x.[__getitem__](#)(y) <=> x[y]

Kwargs:
key (str, int, slice, list): Item to retrieve.

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

__init__(self, recordsInst, recordId, dataType= '/')
Initiate the class.

Kwargs:
recordsInst ([biomesh.records](#)): Instance of the [records](#) collection.
recordId (int): ID of the record.
dataType (str): Type of the data. Default: '/'.

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

__iter__(self)
Iterator operator.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

__len__(self)
Number of datasets in the record.

x.[__len__](#)() <=> len(x)

Kwargs:

Kwrvals:

See Also:

Notes:
If the case (signal or events) is still unspecified, returns the total number of datasets. Otherwise, returns the number of datasets.

Example:

References:
.. [1]

__str__(self)
str operator.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

list(self)
List the datasets.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

class experiments

Class to operate on the Experiments collection.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

Methods defined here:

__init__(self, db=None, meta=None, parent=None, **kwargs)
Initialize the [experiments](#) class.

Kwargs:
db (pymongo.database.Database): An instance representing the DB.
meta (str): The path to the meta HDF5 file.

Kwrvals:

See Also:

Notes:
For [experiments](#), the key 'name' is the preferred indexer.

Example:

References:
.. [1]

add(self, experiment=None, flag=True)
To add an experiment to the DB's '[experiments](#)' collection.

Kwargs:
experiment (dict): Experiment (JSON) to add.
flag (bool): Flag to store in HDF5. Default: True.

Kwrvals:
experimentId (int): ID of the experiment in the DB.

See Also:

Notes:
If the experiment already exists (i.e. there is an experiment with the same name) the ID of the experiment in the DB is returned.

Example:
expId = db.[experiments.add](#)({'name': 'experiment'})['experimentId']

References:
.. [1]

get(self, refine={}, restrict={})
Make a general query the '[experiments](#)' collection.

Kwargs:
refine (dict): Dictionary to refine the search. Default: {}.
restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:
doclist (list): The list of documents with the results of the query. Empty list if no match is found.

See Also:

Notes:

Example:
get all [experiments](#)
res = db.[experiments.get](#)()['docList']
get [experiments](#) with 'field' set to 'new'
res = db.[experiments.get](#)(refine={'field': 'new'})['docList']
get [experiments](#) with 'field' set to 'new' and 'flag' set to False
res = db.[experiments.get](#)(refine={'field': 'new', 'flag': False})['docList']

References:
.. [1]

getById(self, experimentId=None, restrict={})
Query the '[experiments](#)' collection by ID.

Kwargs:
experimentId (int): ID of the experiment to query.
restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:
doc (dict): The document with the results of the query. None if no match is found.

See Also:

Notes:

Example:

```
# get everything
doc = db.experiments.getById(0)['doc']
# only return the ID
doc = db.experiments.getById(0, restrict={'_id': 1})['doc']
# only return the name (ID is also returned!)
doc = db.experiments.getById(0, restrict={'name': 1})['doc']
# only return the name and don't return the ID
doc = db.experiments.getById(0, restrict={'name': 1, '_id': 0})['doc']
```

References:

.. [1]

getByName(self, experimentName=None, restrict={})

Query the '[experiments](#)' collection by name.

Kwargs:

experimentName (str): Name of the experiment to query.

restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:

doc (dict): The document with the results of the query. None if no match in found.

See Also:

Notes:

Example:

```
# get everything
doc = db.experiments.getByName('experiment')['doc']
# only return the ID
doc = db.experiments.getByName('experiment', restrict={'_id': 1})['doc']
# only return the name (ID is also returned!)
doc = db.experiments.getByName('experiment', restrict={'name': 1})['doc']
# only return the name and don't return the ID
doc = db.experiments.getByName('experiment', restrict={'name': 1, '_id': 0})['doc']
```

References:

.. [1]

list(self)

List all the [experiments](#) in the database.

Kwargs:

Kwrvals:

nameList (list): List with the names of the [experiments](#) in the database.

See Also:

Notes:

Example:

```
explist = db.experiments.list()['nameList']
```

References:

.. [1]

update(self, experimentName=None, info={})

Update an experiment with the given information. Fields can be added, and its type changed, but not deleted.

Kwargs:

experimentName (str): Name of the experiment to update.

info (dict): Dictionary with the information to update. Default: {}.

Kwrvals:

See Also:

Notes:

Example:

```
db.experiments.update('experiment', {'new': 'field'})
```

References:

.. [1]

class records

Class to operate on the Records collection.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

Methods defined here:

__getitem__(self, key)

To get Records from the DB.

x.[__getitem__](#)(y) <=> x[y]

Kwargs:

key (int, slice, list): Items to retrieve.

Kwrvals:

See Also:

Notes:

Example:

```
container = db.records[0]
containerList = db.records[:]
containerList = db.records[[0, 3, 4]]
```

References:

.. [1]

__init__(self, db=None, path=None, srvPath=None, meta=None, parent=None, **kwargs)

Initialize the [records](#) class.

Kwargs:

db (pymongo.database.Database): An instance representing the DB.

path (str): The path to store the HDF5 files.

srvPath (str): The path to store the HDF5 files on the remote server.

meta (str): The path to the meta HDF5 file.

Kwrvals:

See Also:

Notes:

For [records](#), the key '_id' is the preferred indexer.

Example:

References:

.. [1]

__iter__(self)

Iterator operator over the Records in the DB.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

```
for item in db.records:
    signal = item['signals']['test'][0].signal
    metadataS = item['signals']['test'][0].metadata
    timeStamps = item['events']['test'][0].timeStamps
    values = item['events']['test'][0].values
    metadataE = item['events']['test'][0].metadata
```

References:

.. [1]

__len__(self)

Returns the number of [records](#) on the DB.

x.[__len__](#)() <=> len(x)

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

len(db.[records](#))

References:

.. [1]

add(self, record=None, flag=True)

To add a record to the DB's '[records](#)' collection.

Kwargs:

record (dict): Record (JSON) to add.

flag (bool): Flag to store in HDF5. Default: True.

Kwrvals:

recordId (int): ID of the record in the DB.

See Also:

Notes:

A new record is always added, regardless of the existence of [records](#) with the same name. If the record to add has no 'experiment' ar

Example:

recId = db.[records.add](#)({'name': 'record'})['recordId']

References:

.. [1]

addAudit(self)

To add information to 'audit' field.

TO DO!

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

addEvent(self, recordId=None, eventType='/', timeStamps=None, values=None, mdata={}, flag=True, compress=False)

To add events (asynchronous data) to a record.

Kwargs:

recordId (int): ID of the record.

eventType (str): Type of the events to add. Default: '/'.

timeStamps (array): Array of time stamps. Default: [].

values (array): Array with data for each time stamp. Default: [].

mdata (dict): JSON with metadata about the events. Default: {}.

flag (bool): Flag to store in HDF5. Default: True.

compress (bool): Flag to compress the data (GZIP). Default: False.

Kwrvals:

recordId (int): ID of the record.

eventRef (str): Storage name of the events.

eventType (str): Type of the inserted events.

See Also:

Notes:

Example:

```
res = db.records.addEvent(0, '/test', [0, 1, 2], [[0, 1], [2, 3], [4, 5]], {'comments': 'test event'})
```

References:

.. [1]

addSignal(self, recordId=None, signalType='/', signal=None, mdata={}, flag=True, compress=False, updateDuration=False)
To add signal (synchronous) data to a record.

Kwargs:

recordId (int): ID of the record.

signalType (str): Type of the signal to add. Default: '/'.

signal (array): Signal to add. Default: [].

mdata (dict): Dictionary (JSON) with metadata about the signal. Default: {}.

flag (bool): Flag to store in HDF5. Default: True.

compress (bool): Flag to compress the data (GZIP). Default: False.

updateDuration (bool): Flag to update the duration of the record. Default: False.

Kwrvals:

recordId (int): ID of the record.

signalRef (str): Storage name of the signal.

signalType (str): Type of the inserted signal.

See Also:

Notes:

Example:

```
res = db.records.addSignal(0, '/test', [0, 1, 2, 3], {'comments': 'test signal'})
```

References:

.. [1]

addTags(self, recordId=None, tags=None, flag=True)
Add tags to a record.

Kwargs:

recordId (int): ID of the record.

tags (list): Tags to add.

flag (bool): Flag to store in HDF5. Default: True.

Kwrvals:

None

See Also:

Notes:

Example:

```
db.records.addTags(0, ['a', 'b', 'c'])
```

References:

.. [1]

delEvent(self, recordId=None, eventType='/', eventRef=None)
Remove events from a record.

Kwargs:

recordId (int): ID of the record.

eventType (str): Type of the desired event. Default: '/'.

eventRef (str, int): Storage name (global) or index (local) of the desired events.

Kwrvals:

See Also:

Notes:

Example:

```
db.records.delEvent(0, '/test', 0)
```

References:

.. [1]

delEventType(self, recordId=None, eventType='/')

Remove an event type from a record (including all sub-types).

Kwargs:

recordId (int): ID of the record.

eventType (str): Type of the desired event. Default: '/'.

Kwrvals:

See Also:

Notes:

Example:

db.[records.delEventType](#)(0, '/test')

References:

.. [1]

delSignal(self, recordId=None, signalType='/', signalRef=None)

Remove a signal (synchronous data) from a record.

Kwargs:

recordId (int): ID of the record.

signalType (str): Type of the desired data. Default: '/'.

signalRef (str, int): Storage name (global) or index (local) of the desired signal.

Kwrvals:

See Also:

Notes:

Example:

db.[records.delSignal](#)(0, '/test', 0)

References:

.. [1]

delSignalType(self, recordId=None, signalType='/')

Remove a signal type from a record (including all sub-types).

Kwargs:

recordId (int): ID of the record.

signalType (str): Type of the desired data. Default: '/'

Kwrvals:

See Also:

Notes:

Example:

db.[records.delSignalType](#)(0, '/test')

References:

.. [1]

delTags(self, recordId=None, tags=None)

Delete tags from a record.

Kwargs:

recordId (int): ID of the record.

tags (list): Tags to add.

Kwrvals:

None

See Also:

Notes:

Example:

db.[records.delTags](#)(0, ['a', 'b', 'c'])

References:

.. [1]

delete(self, recordId=None, keepFile=True)

Remove a record from the database.

Kwargs:

recordId (int): ID of the record.

keepFile (bool): Flag to keep local file. Default: True.

Kwrvals:

See Also:

Notes:

Example:

```
db.records.delete(0)
```

References:

.. [1]

get(self, refine={}, restrict={})

Make a general query the '[records](#)' collection.

Kwargs:

refine (dict): Dictionary to refine the search. Default: {}.

restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:

doclist (list): The list of documents with the results of the query. Empty list if no match is found.

See Also:

Notes:

Example:

```
# get all records
res = db.records.get()['docList']
# get records with 'field' set to 'new'
res = db.records.get(refine={'field': 'new'})['docList']
# get records with 'field' set to 'new' and 'flag' set to False
res = db.records.get(refine={'field': 'new', 'flag': False})['docList']
```

References:

.. [1]

getAll(self, restrict={}, count=-1, randomFlag=False)

Generate a list of [records](#) present in the DB. The list may include all [records](#), [records](#) pertaining to an experiment (or list of [experiments](#))

Kwargs:

restrict (dict): To restrict the results. Can have the following keys:

experiment (int, str, list): Experiment ID, experiment name, list of experiment IDs, list of experiment names, or list of experiments

subject (int, str, list): Subject ID, subject name, list of subject IDs, list of subject names, or list of experiment names and subjects

Default: {}

count (int): The resulting list has, at most, 'count' items. Set to -1 to output all [records](#) found. Default: -1.

randomFlag (bool): Set this flag to True to randomize the output list. Default: False.

Kwrvals:

idList (list): List with the [records](#) IDs that match the search.

See Also:

Notes:

Example:

```
# get all records
idList = db.records.getAll()['idList']
# get at most 5 records
idList = db.records.getAll(count=5)['idList']
# get all records from experiment 'exp'
idList = db.records.getAll(restrict={'experiment': 'exp'})['idList']
# get all records from experiment 'exp' and 'dexp'
idList = db.records.getAll(restrict={'experiment': ['exp', 'dexp']})['idList']
# get all records from subject 0
idList = db.records.getAll(restrict={'subject': 0})['idList']
# get all records from subject 0 and 1
idList = db.records.getAll(restrict={'subject': [0, 1]})['idList']
# get all records from subject 0 and 1 and 'John Smith'
idList = db.records.getAll(restrict={'subject': [0, 1, 'John Smith']})['idList']
# get all records from experiment 'exp', but only from subjects 0 and 1
idList = db.records.getAll(restrict={'experiment': 'exp', 'subject': [0, 1]})['idList']
```

References:

.. [1]

getId(self, recordId=None, restrict={})

Query the '[records](#)' collection by ID.

Kwargs:
recordId (int): ID of the record.

restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:
doc (dict): The document with the results of the query. None if no match in found.

See Also:

Notes:

Example:

```
# get everything
doc = db.records.getById(0)['doc']
# only return the ID
doc = db.records.getById(0, restrict={'_id': 1})['doc']
# only return the name (ID is also returned!)
doc = db.records.getById(0, restrict={'name': 1})['doc']
# only return the name and don't return the ID
doc = db.records.getById(0, restrict={'name': 1, '_id': 0})['doc']
```

References:
.. [1]

getName(self, recordName=None, restrict={})
Query the 'records' collection by name.

Kwargs:
recordName (str): Name of the record.

restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:
doc (dict): The document with the results of the query. None if no match in found.

See Also:

Notes:

Example:

```
# get everything
doc = db.records.getName('record')['doc']
# only return the ID
doc = db.records.getName('record', restrict={'_id': 1})['doc']
# only return the name (ID is also returned!)
doc = db.records.getName('record', restrict={'name': 1})['doc']
# only return the name and don't return the ID
doc = db.records.getName('record', restrict={'name': 1, '_id': 0})['doc']
```

References:
.. [1]

getEvent(self, recordId=None, eventType='/', eventRef=None)
Retrieve events (asynchronous data) from a record.

Kwargs:
recordId (int): ID of the record.

eventType (str): Type of the desired event. Default: '/'.

eventRef (str, int): Storage name (global) or index (local) of the desired events.

Kwrvals:
timeStamps (array): Array of time stamps.

values (array): Array with data for each time stamp.

mdata (dict): Dictionary with metadata about the events.

See Also:

Notes:

Example:

```
res = db.records.getEvent(0, '/test', 0)
timeStamps = res['timeStamps']
values = res['values']
metadata = res['mdata']
```

References:
.. [1]

getSignal(self, recordId=None, signalType='/', signalRef=None)
Retrieve a signal(synchronous data) from a record and corresponding metadata.

Kwargs:
recordId (int): ID of the record.

signalType (str): Type of the desired signal. Default: '/'.

signalRef (str, int): Storage name (global) or index (local) of the desired signal.

signalRef (int, str): Storage name (global) or index (local) of the desired signal.

Kwrvals:

- signal (array): Array with the signal.
- mdata (dict): Dictionary with the signal's accompanying metadata.

See Also:

Notes:

Example:

```
res = db.records.getSignal(0, '/test', 0)
signal = res['signal']
metadata = res['mdata']
```

References:

.. [1]

listAndTags(self, tags=None)

Lists the [records](#) that simultaneously have all the given tags (AND operator).

Kwargs:

- tags (list): Tags to match.

Kwrvals:

- idList (list): List with the [records](#)' IDs that match the search.

See Also:

Notes:

Example:

```
res = db.records.listAndTags(['a', 'b', 'c'])['idList']
```

References:

.. [1]

listNotTags(self, tags=None)

listOrTags(self, tags=None)

Lists the [records](#) that have, at least, one of the given tags (OR operator).

Kwargs:

- tags (list): Tags to match.

Kwrvals:

- idList (list): List with the [records](#)' IDs that match the search.

See Also:

Notes:

Example:

```
res = db.records.listOrTags(['a', 'b', 'c'])['idList']
```

References:

.. [1]

listSymbolicTags(self, query=None)

listTypes(self, recordId=None)

To list the types of signals and events.

Kwargs:

- recordId (int): ID of the record.

Kwrvals:

- signalTypes (list): List of data types.
- eventTypes (list): List of event types.

See Also:

Notes:

Example:

```
res = db.records.listTypes(0)
```

References:

.. [1]

update(self, recordId=None, info={}, flag=True)

Update a record with the given information. Fields can be added, and its type changed, but not deleted.

Kwargs:

- recordId (int): ID of the record to update.

info (dict): Dictionary with the information to update. Default: {}.

flag (bool): Flag to store in HDF5. Default: True.

Kwrvals:

See Also:

Notes:

Example:

```
db.records.update(0, {'new': 'field'})
```

References:

.. [1]

class **subjects**

Class to operate on the Subjects collection.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

Methods defined here:

__init__(self, db=None, meta=None, parent=None, **kwargs)
Initialize the [subjects](#) class.

Kwargs:

db (pymongo.database.Database): An instance representing the DB.

meta (str): The path to the meta HDF5 file.

Kwrvals:

See Also:

Notes:

For [subjects](#), the key '_id' is the preferred indexer.

Example:

References:

.. [1]

add(self, subject=None, flag=True)
To add a subject to the '[subjects](#)' collection.

Kwargs:

subject (dict): Subject (JSON) to add.

flag (bool): Flag to store in HDF5. Default: True.

Kwrvals:

subjectId (int): ID of the subject in the DB.

See Also:

Notes:

If the subject already exists (i.e. there is a subject with the same name) the ID of the subject in the DB is returned.

Example:

```
subId = db.subjects.add({'name': 'subject'})['subjectId']
```

References:

.. [1]

get(self, refine={}, restrict={})
Make a general query to the '[subjects](#)' collection.

Kwargs:

refine (dict): Dictionary to refine the search. Default: {}.

restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:

doclist (list): The list of dictionaries with the results of the query. Empty list if no match is found.

See Also:

Notes:

Example:

```
# get all subjects
res = db.subjects.get()['doclist']
# get subjects with 'field' set to 'new'
res = db.subjects.get(refine={'field': 'new'})['doclist']
# get subjects with 'field' set to 'new' and 'flag' set to False
res = db.subjects.get(refine={'field': 'new', 'flag': False})['doclist']
```

References:

.. [1]

getById(self, subjectId=None, restrict={})

Query the '[subjects](#)' collection by ID.

Kwargs:

subjectId (int): ID of the subject to query.

restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:

doc (dict): The document with the results of the query. None if no match is found.

See Also:

Notes:

Example:

```
# get everything
doc = db.subjects.getById(0)['doc']
# only return the ID
doc = db.subjects.getById(0, restrict={'_id': 1})['doc']
# only return the name (ID is also returned!)
doc = db.subjects.getById(0, restrict={'name': 1})['doc']
# only return the name and don't return the ID
doc = db.subjects.getById(0, restrict={'name': 1, '_id': 0})['doc']
```

References:

.. [1]

getByName(self, subjectName=None, restrict={})

Query the '[subjects](#)' collection by name.

Kwargs:

subjectName (str): Name of the subject to query.

restrict (dict): Dictionary to restrict the information sent by the DB. Default: {}.

Kwrvals:

doc (dict): The document with the results of the query. None if no match is found.

See Also:

Notes:

Example:

```
# get everything
doc = db.subjects.getByName('subject')['doc']
# only return the ID
doc = db.subjects.getByName('subject', restrict={'_id': 1})['doc']
# only return the name (ID is also returned!)
doc = db.subjects.getByName('subject', restrict={'name': 1})['doc']
# only return the name and don't return the ID
doc = db.subjects.getByName('subject', restrict={'name': 1, '_id': 0})['doc']
```

References:

.. [1]

list(self)

List all the [subjects](#) in the database.

Kwargs:

Kwrvals:

idList (list): List with the IDs of the [subjects](#) in the database.

See Also:

Notes:

Example:
subList = db.[subjects.list](#)()['idList']

References:
.. [1]

update(self, subjectId=None, info={})

Update a subject with the given information. Fields can be added, and its type changed, but not deleted.

Kwargs:
subjectId (int): ID of the subject to update.

info (dict): Dictionary with the information to update. Default: {}.

Kwrvals:

See Also:

Notes:

Example:
db.[subjects.update](#)(0, {'new': 'field'})

References:
.. [1]

Functions

listDB(host='localhost', port=27017)

Lists the databases present in a given server.

Kwargs:
host (str): Network address of the MongoDB server. Default: 'localhost'.

port (int): Port the MongoDB server is listening on. Default: 27017.

Kwrvals:
dbList (list): A list with DB names.

See Also:

Notes:

Example:
[listDB](#)('193.136.222.234', 27017)

References:
.. [1]

Data

RETAG = <_sre.SRE_Pattern object>

```
.. module:: philipsXML
   :platform: Windows
   :synopsis: This module provides tools to read ECG records stored using the Philips SierraECG standard.

.. moduleauthor:: Carlos Carreiras
```

Modules

[lxml.etree](#) [numpy](#) [xmldict](#)
[bspyclone.database.h5db](#) [os](#)

Functions

```
readXML(fpath=None, validate=True)
    Read a decompressed Philips XML file.

    Kwargs:
        fpath (str): XML file to read.

        validate (bool): If True, validates the XML schema.

    Kwrvals:
        date (str): ISO 8601 acquisition date and time.

        labels (list): Label for each acquired ECG channel.

        resolution (int): Signal resolution.

        sampleRate (float): Acquisition sample rate.

        signal (array): Signal array, where each line is an ECG channel.

        subjectID (str): The ID of the subject.

        subjectName (str): The name of the subject.

        xmlFileName (str): Name of the source XML file.

    See Also:

    Notes:

    Example:

    References:
        .. [1]
```

[bspyclone.database.physionet](#) [index](#) [d:\work\productioncode\clones\biosppy\bspyclone\database\physionet.py](#)

```
.. module:: biomesb
    :platform: Unix, Windows
    :synopsis: Convert Physionet database files to csv (requires wfdb).

.. moduleauthor:: Carlos Carreiras
```

Modules

csv	numpy	subprocess
fnmatch	os	

Functions

```
batchConvert(basePath, db, dirPath=None)

convert2csv(basePath, record, dirPath)

convertAnnotation2CSV(basePath, record, dirPath)

listFiles(basePath, db)

readCSV(path)
```

```
.. module:: ecg
    :platform: Unix, Windows
    :synopsis: This module provides various functions to handle ECG signals.

.. moduleauthor:: Filipe Canento
```

Modules

[bspyclone.ecg.models](#)
[bspyclone.filt](#)

[bspyclone.plux](#)
[pylab](#)

[scipy](#)
[bspyclone.tools](#)

[unittest](#)

Functions

ecg(Signal=None, SamplingRate=1000.0, Filter={})
Determine ECG signal information.

Kwargs:

Signal (array): input ECG signal.

SamplingRate (float): Sampling frequency (Hz).

Filter (dict): Filter parameters.

Kwrvals:

R (array): heart beat indexes (or instants in seconds if sampling rate is defined).

Configurable fields:{"name": "ecg.ecg", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter"], "outputs": ["R"]}

See Also:

filt
models.hamilton

Notes:

Example:

References:

.. [1]

features(Signal=None, SamplingRate=1000.0, Filter={})
Retrieves relevant ECG signal features.

Kwargs:

Signal (array): input ECG signal.

SamplingRate (float): Sampling frequency (Hz).

Filter (dict): Filter parameters.

Kwrvals:

R (array): ECG R-peak indexes (or instants in seconds if sampling rate is defined)

mean (float): mean

std (float): standard deviation

var (float): variance

skew (ndarray): skewness

kurtosis (array): kurtosis

ad (float): absolute deviation

Configurable fields:

{"name": "ecg.features", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter"], "outputs": ["R", "mean", "std", "var", "skew",

See Also:

filt
models.hamilton
tls.statsf

Notes:

Example:

References:

.. [1]

filt(Signal=None, SamplingRate=1000.0, UpperCutoff=16.0, LowerCutoff=8.0, Order=4)
Filters an input ECG signal.

By default, the return is the filtered Signal assuming a

1000Hz sampling frequency and the following filter sequence:

1. 4th order low-pass filter with cutoff frequency of 16Hz;
2. 4th order high-pass filter with cutoff frequency of 8Hz;
3. d[]/dt;

4. 80ms Hamming Window Smooth.

Kwargs:

Signal (array): input signal.

SamplingRate (float): Sampling frequency (Hz).

UpperCutoff (float): Low-pass filter cutoff frequency (Hz).

LowerCutoff (float): High-pass filter cutoff frequency (Hz).

Order (int): Filter order.

Kwrvals:

Signal (array): output filtered signal.

Configurable fields:

{"name": "ecg.filt", "config": {"UpperCutoff": "16.", "SamplingRate": "1000.", "LowerCutoff": "8.", "Order": "4"}, "inputs": ["Signal"], "out

See Also:

`flt.zpdf`

`flt.smooth`

Notes:

Example:

```
Signal = load(...)
```

```
SamplingRate = ...
```

```
res = filt(Signal=Signal, SamplingRate=SamplingRate)
```

```
plot(res['Signal'])
```

References:

.. [1] P.S. Hamilton, Open Source ECG Analysis Software Documentation, E.P.Limited
<http://www.eplimited.com/osea13.pdf>


```
.. module:: hrv
   :platform: Unix, Windows
   :synopsis: This module provides various methods to perform Heart Rate Variability analysis.

.. moduleauthor:: Carlos Carreiras
```

Modules

[numpy](#)[scipy.signal](#)

Functions

hrv(R, sampleRate, npoints=100000)

Compute Heart Rate Variability (HRV) metrics from a sequence of R peak locations.

Kwargs:

R (list, array): Positions of the R peaks in samples (from a segmentation algorithm).

sampleRate (float): Sampling rate (Hz).

npoints (int): Number of frequency points for the Lomb-Scargle periodogram. Default=100000.

Kwrvals:

time (array): Instantaneous heart rate time points.

HR (array): Instantaneous heart rate.

RR (array): Instantaneous RR intervals.

mHR (float): Mean heart rate.

SD (float): Standard deviation of heart rate.

RMSSD (float): Root mean square of successive differences.

HF (float): High frequency power (0.15 to 0.4 Hz).

LF (float): Low frequency power (0.04 to 0.15 Hz).

L2HF (float): Ratio of LF to HF power.

See Also:

Notes:

Maybe deal with un/mis-detected R peaks

Example:

References:

.. [1]

```
.. module:: models
    :platform: Unix, Windows
    :synopsis: This module provides various methods to segment ECG signals.

.. moduleauthor:: Filipe Canento, Carlos Carreiras, Francisco David
```

Modules

collections	pylab	scipy.signal	traceback
numpy	scipy	bspyclone.ecg.tools	

Functions

ESSF(Signal=None, SamplingRate=1000.0, ESSF_params={'s_amp': 70.0, 's_win': 0.3})
 ECG Slope Sum Function: algorithm to detect ECG beat indexes.

Kwargs:

- Signal (array): input filtered ECG signal.
- SamplingRate (float): Sampling frequency (Hz).
- ESSF_params (float): data dependent parameters.

Returns:

- rpeaks (array): R-peak indexes

Configurable fields:{"name": "models.ESSF", "config": {"SamplingRate": "1000."}, "inputs": ["Signal"], "outputs": ["rpeaks"]}

See Also:

Notes:

Example:

References:

- .. [1] ...

armSSF(Signal=None, SamplingRate=1000.0, threshold=20, winB=0.03, winA=0.01, Params=None)
 Modified Slope Sum Function R peak detection algorithm for ARM board.

Kwargs:

- Signal (array): input ECG signal.
- SamplingRate (float): Sampling frequency (Hz).
- threshold (int): threshold.
- winB (int): size of search window before candidate.
- winA (int): size of search window after candidate.

Returns:

- rpeaks (array): R-peak indexes.

Configurable fields: {}

See Also:

Notes:

Example:

References:

- .. [1] ...

batch_christov(Signal=None, SamplingRate=1000.0, debug=False, IF=True)

batch_engzee(Signal=None, SamplingRate=1000.0, debug=False, IF=True)

Kwargs:

Kwrvals:

Configurable fields:{"name": "models.batch_engzee", "config": {"SamplingRate": "1000.0"}, "inputs": ["Signal"], "outputs": []}

See Also:

Notes:

Example:

References:

- .. [1]

christov(Signal=None, SamplingRate=1000.0, Filter={}, Params={}, Show=False)
Determine ECG signal information.

Kwargs:

Signal (array): input ECG signal.

SamplingRate (float): Sampling frequency (Hz).

Filter (dict): Filter parameters.

Params (dict): Initial conditions.

Kwrvals:

Configurable fields:{"name": "models.christov", "config": {"SamplingRate": "1000.", "Show": "False"}, "inputs": ["Signal", "Filter", "Params"]

See Also:

Notes:

Example:

References:

- .. [1] Ivaylo I Christov, Real time electrocardiogram QRS detection using combined adaptive threshold, BioMedical Engineering OnLine 2004, 3:28
This article is available from: <http://www.biomedical-engineering-online.com/content/3/1/28>
2004 Christov; licensee BioMed Central Ltd.

definepeak(signal, srate)

engzee(Signal=None, SamplingRate=1000.0, Params=None)

Determine ECG signal information. Adaptation of the Engelse and Zeelenberg by [1].

Kwargs:

Signal (array): input ECG signal.

SamplingRate (float): Sampling frequency (Hz).

Params (dict): Initial conditions:

CSignal (array): Continuity signal

MM (array): MM threshold

MMidx (int): MM threshold current index

NN (array): NN threshold

NNidx (int): NN threshold current index

offset (array): last intersection points

prevR (collections.deque): Last R positions

Rminus (float): Segmentation window size to the left of R (default: 200 ms)

Rplus (float): Segmentation window size to the right of R (default: 400 ms)

update (bool): Update flag

Kwrvals:

R (array): R peak indexes

Segments (array): Extracted segments

HH (float): Heart rate

Params (dict):

CSignal (array): Continuity signal

MM (array): MM threshold

MMidx (int): MM threshold current index

NN (array): NN threshold

NNidx (int): NN threshold current index

offset (array): last intersection points

prevR (collections.deque): Last R positions

Rminus (float): Segmentation window size to the left of R (default: 200 ms)

Rplus (float): Segmentation window size to the right of R (default: 400 ms)

update (bool): Update flag

See Also:

Notes:

Example:

References:

```
.. [1] Andre Lourenco, Hugo Silva, Paulo Leite, Renato Lourenco and Ana Fred,  
REAL TIME ELECTROCARDIOGRAM SEGMENTATION FOR FINGER BASED ECG BIOMETRICS
```

engzee_incomplete(Signal=None, SamplingRate=1000.0, Filter=None, Params=None, initialFilter=False)

--- REVIEW DOCS ---

Determine ECG signal information. Adaptation of the Engelse and Zeelenberg by [1].

Kwargs:

Signal (array): input ECG signal.

SamplingRate (float): Sampling frequency (Hz).

Filter (dict): Filter parameters.

Params (dict): Initial conditions

Segments (list): list of ECG segments

RawSegments (list): list of Raw ECG segments

MM (array): MM threshold

MMidx (int): MM threshold current index

offset (array): last intersection points

offsetidx (int): offset current index

rpeak (array): R indexes

nthfpluss (array): nthfplus intersection points

Rminus (int): Segmentation window size to the left of R (default: 200 ms)

Rplus (int): Segmentation window size to the right of R (default: 400 ms)

initialFilter (bool): Apply hands ECG filter (default: False).

Kwrvals:

R (array): R peak indexes

Params (dict):

Segments (list): list of ECG segments

RawSegments (list): list of Raw ECG segments

MM (array): MM threshold

MMidx (int): MM threshold current index

offset (array): last intersection points

offsetidx (int): offset current index

rpeak (array): R indexes

nthfpluss (array): nthfplus intersection points

Rminus (int): Segmentation window size to the left of R (default: 200 ms)

Rplus (int): Segmentation window size to the right of R (default: 400 ms)

Configurable fields:

{**"name"**: "models.engzee", **"config"**: {"SamplingRate": "1000.", "initialfilter": "True"}, **"inputs"**: ["Signal", "Filter", "Params"], **"outputs"**:

See Also:

Notes:

1. Tested for chest ECGs: paper thresholds, no initial filtering
2. Finger/palm ECGs: there's an initial filtering and some parameters were modified

Example:

Data filename

fname = ...

Load data

dataecg = plux.loadbpf(fname, usecols=(3,))

Init

mrkr = 0

SamplingRate = dataecg.header['SamplingFrequency']

win = 3*SamplingRate

mrkrend = win

ECG algorithm, segmentation: R-200ms to R+400ms

```

res = ecgmodule.models.engzee(Signal=dataecg[:win], SamplingRate=SamplingRate, Params={'Rminus': int(0.2*SamplingRate), 'Rplus': int(0.4*
dlen = len(dataecg)
while ( mrkrend < dlen ):
    mrkr = res['Params']['offset'][-1]
    res = ecgmodule.models.engzee(Signal=dataecg[mrkr:mrkrend], SamplingRate=SamplingRate, Params=res['Params'])
ECGSegments = res['Params']['Segments']

```

References:

.. [1] Andre Lourenco, Hugo Silva, Paulo Leite, Renato Lourenco and Ana Fred,
REAL TIME ELECTROCARDIOGRAM SEGMENTATION FOR FINGER BASED ECG BIOMETRICS

engzee_old(Signal=None, SamplingRate=1000.0, initialfilter=True, Filter={}, Params={})
HAS PROBLEM WITH R PEAKS LOCATION IN ONLINE USE-CASE

Determine ECG signal information. Adaptation of the Engelse and Zeelenberg by [1].

Kwargs:

Signal (array): input ECG signal.

SamplingRate (float): Sampling frequency (Hz).

Filter (dict): Filter parameters.

Params (dict): Initial conditions

Segments (list): list of ECG segments

RawSegments (list): list of Raw ECG segments

MM (array): MM threshold

MMidx (int): MM threshold current index

offset (array): last intersection points

offsetidx (int): offset current index

rpeak (array): R indexes

nthfpluss (array): nthfplus intersection points

Rminus (int): Segmentation window size to the left of R (default: 200 ms)

Rplus (int): Segmentation window size to the right of R (default: 400 ms)

Kwrvals:

R (array): R peak indexes

Params (dict):

Segments (list): list of ECG segments

RawSegments (list): list of Raw ECG segments

MM (array): MM threshold

MMidx (int): MM threshold current index

offset (array): last intersection points

offsetidx (int): offset current index

rpeak (array): R indexes

nthfpluss (array): nthfplus intersection points

Rminus (int): Segmentation window size to the left of R (default: 200 ms)

Rplus (int): Segmentation window size to the right of R (default: 400 ms)

Configurable fields:

{"name": "models.engzee", "config": {"SamplingRate": "1000.", "initialfilter": "True"}, "inputs": ["Signal", "Filter", "Params"], "outputs":

See Also:

Notes:

1. Tested for chest ECGs: paper thresholds, no initial filtering
2. Finger/palm ECGs: there's an initial filtering and some parameters were modified

Example:

```

# Data filename
fname = ...

# Load data
dataecg = plux.loadbpf(fname, usecols=(3,))

# Init

```

```

mrkr = 0

SamplingRate = dataecg.header['SamplingFrequency']

win = 3*SamplingRate

mrkrend = win

# ECG algorithm, segmentation: R-200ms to R+400ms

res = ecgmodule.models.engzee(Signal=dataecg[:win], SamplingRate=SamplingRate, Params={'Rminus': int(0.2*SamplingRate), 'Rplus': int(0.4*SamplingRate)})

dlen = len(dataecg)

while ( mrkrend < dlen ):

    mrkr = res['Params']['offset'][-1]

    res = ecgmodule.models.engzee(Signal=dataecg[mrkr:mrkrend], SamplingRate=SamplingRate, Params=res['Params'])

    mrkrend = res['Params']['Rplus']

ECGSegments = res['Params']['Segments']

```

References:

```

.. [1] Andre Lourenco, Hugo Silva, Paulo Leite, Renato Lourenco and Ana Fred,
REAL TIME ELECTROCARDIOGRAM SEGMENTATION FOR FINGER BASED ECG BIOMETRICS

```

gamboa(Signal=None, SamplingRate=1000.0, tol=0.002)

Gamboa's algorithm to detect ECG beat indexes.

Kwargs:

```

Signal (array): input filtered ECG signal.

SamplingRate (float): Sampling frequency (Hz).

tol (float): tolerance parameter.

```

Returns:

```

rpeaks (array): R-peak indexes

```

Configurable fields:{"name": "models.gamboa", "config": {"SamplingRate": "1000."}, "inputs": ["Signal"], "outputs": ["rpeaks"]}

See Also:

Notes:

Example:

References:

```

.. [1] ...

```

hamilton(Signal=None, SamplingRate=1000.0, Params=None, hand=True)

Algorithm to detect ECG beat indexes.

Kwargs:

```

Signal (array): input filtered ECG signal.

SamplingRate (float): Sampling frequency (Hz).

Params (dict): Initial conditions:

hand (bool): Signal is obtained from the hands.

```

Kwrvals:

```

Signal (array): output filtered signal if Filter is defined.

R (array): R peak indexes (or instants in seconds if sampling rate is defined).

init (dict): dict with initial values of some variables

    npeaks (int): number of detected heart beats.

    indexqrs (int): most recent QRS complex index.

    indexnoise (int): most recent noise peak index.

    indexrr (int): most recent R-to-R interval index.

    qrspeakbuffer (array): 8 most recent QRS complexes.

    noisepeakbuffer (array): 8 most recent noise peaks.

    rrinterval (array): 8 most recent R-to-R intervals.

    DT (float): QRS complex detection threshold.

    offset (int): signal start in samples.

```

Configurable fields:

```

{"name": "models.hamilton", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter", "init"], "outputs": ["Signal", "R", "init", "DT"]}

```

See Also:

```

filt

```

Notes:

Example:

References:

.. [1] P.S. Hamilton, Open Source ECG Analysis Software Documentation, E.P.Limited
<http://www.eplimited.com/osea13.pdf>

hamilton_old(Signal=None, SamplingRate=1000.0, Filter=(), init=())
Algorithm to detect ECG beat indexes.

Kwargs:

Signal (array): input filtered ECG signal.
SamplingRate (float): Sampling frequency (Hz).
Filter (dict): Filter parameters.

Kwrvals:

Signal (array): output filtered signal if Filter is defined.
R (array): R peak indexes (or instants in seconds if sampling rate is defined).
init (dict): dict with initial values of some variables
 npeaks (int): number of detected heart beats.
 indexqrs (int): most recent QRS complex index.
 indexnoise (int): most recent noise peak index.
 indexrr (int): most recent R-to-R interval index.
 qrspeakbuffer (array): 8 most recent QRS complexes.
 noisepeakbuffer (array): 8 most recent noise peaks.
 rrinterval (array): 8 most recent R-to-R intervals.
 DT (float): QRS complex detection threshold.
 offset (int): signal start in samples.

Configurable fields:

{"name": "models.hamilton", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter", "init"], "outputs": ["Signal", "R", "init"], '}

See Also:

filt

Notes:

Example:

References:

.. [1] P.S. Hamilton, Open Source ECG Analysis Software Documentation, E.P.Limited
<http://www.eplimited.com/osea13.pdf>

hamilton_tst(Signal=None, SamplingRate=1000.0, init=None, hand=True)
Algorithm to detect ECG beat indexes.

Kwargs:

Signal (array): input filtered ECG signal.
SamplingRate (float): Sampling frequency (Hz).
Filter (dict): Filter parameters.

Kwrvals:

Signal (array): output filtered signal if Filter is defined.
R (array): R peak indexes (or instants in seconds if sampling rate is defined).
init (dict): dict with initial values of some variables
 npeaks (int): number of detected heart beats.
 indexqrs (int): most recent QRS complex index.
 indexnoise (int): most recent noise peak index.
 indexrr (int): most recent R-to-R interval index.
 qrspeakbuffer (array): 8 most recent QRS complexes.
 noisepeakbuffer (array): 8 most recent noise peaks.
 rrinterval (array): 8 most recent R-to-R intervals.
 DT (float): QRS complex detection threshold.
 offset (int): signal start in samples.

Configurable fields:

{"name": "models.hamilton", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter", "init"], "outputs": ["Signal", "R", "init"], '}

See Also:

filt

Notes:

Example:

References:

.. [1] P.S. Hamilton, Open Source ECG Analysis Software Documentation, E.P.Limited
<http://www.eplimited.com/osea13.pdf>

happee(Signal=None, SamplingRate=1000.0, window=3.0, overlap=0.75, threshold=0.0, gridMem=None)
Implementation of the Highly Accurate Pulse Predictor for Exercise Equipment (HAPPEE) [1].

Kwargs:

Signal (array): input ECG signal.

SamplingRate (float): Sampling frequency (Hz).

window (float): Size (seconds) of the search window.

overlap (float): Percentage of window overlap.

threshold (float): Pulse detection threshold.

gridMem (dict): helper grid dictionary.

Returns:

Configurable fields: {}

See Also:

Notes:

Example:

References:

.. [1] Kay, S.; Ding, Q.; Li, D., "On the Design and Implementation of A Highly Accurate Pulse Predictor for Exercise Equipment," IEEE Trans. on Biomedical Engineering, 10.1109/TBME.2015.2407155

happee_grid(length, grid, M)
Generate the helper grid dictionary for the HAPPEE algorithm [1].

Kwargs:

length (int): size of the signal window.

grid (array, list): search grid of pulse periods.

M (int): Pulse length.

Returns:

gridMem (dict): helper grid dictionary.

Configurable fields: {}

See Also:

Notes:

Example:

References:

.. [1] Kay, S.; Ding, Q.; Li, D., "On the Design and Implementation of A Highly Accurate Pulse Predictor for Exercise Equipment," IEEE Trans. on Biomedical Engineering, 10.1109/TBME.2015.2407155

happee_search(window, grid, M, gridMem)
Perform HAPPEE [1] search on a signal window.

Kwargs:

window (array): input signal window.

grid (array, list): search grid of pulse periods.

M (int): Pulse length.

gridMem (dict): helper grid dictionary.

Returns:

T (float): Test statistic.

n0 (int): Location of first pulse.

P (int): Pulse period.

Configurable fields: {}

See Also:

Notes:

Example:

References:

.. [1] Kay, S.; Ding, Q.; Li, D., "On the Design and Implementation of A Highly Accurate Pulse Predictor for Exercise Equipment," IEEE Trans. on Biomedical Engineering, 10.1109/TBME.2015.2407155

happee_variance(signal)

Compute HAPPEE variance [1].

Kwargs:

signal (array): input signal.

Returns:

variance (float): HAPPEE variance.

Configurable fields: {}

See Also:

Notes:

Example:

References:

.. [1] Kay, S.; Ding, Q.; Li, D., "On the Design and Implementation of A Highly Accurate Pulse Predictor for Exercise Equipment," IEEE Trans. on Biomedical Engineering, 10.1109/TBME.2015.2407155

monhe(Signal=None, SamplingRate=1000.0)

GREAT but:

discard crossings close to one another by less than 100 ms

monhe2(Signal=None, Filtered=None, SamplingRate=1000.0, checkSign=True)

GREAT but:

discard crossings close to one another by less than 100 ms

Created on May 28, 2013

@author: Carlos

Modules

[numpy](#)

[pylab](#)

[scipy.signal](#)

Functions

PLFCorrSynchronize(a, v, flag=False)

PLFSynchronize(a, v, flag=False, minOverlap=1)

checkECG(data)

compareSegmentation(referenceR=None, testR=None, SamplingRate=None, offset=0, minRR=None, tol=0.05)

extractHeartbeats(signal, R, sampleRate, before=0.2, after=0.4)

getSyncSlices(d, a, v)

syncOverlap(a, v)

synchronize(a, v, flag=False)

```
.. module:: eda
    :platform: Unix, Windows
    :synopsis: This module provides various functions to handle EDA signals.

.. moduleauthor:: Filipe Canento
```

Modules

[bspyclone.eda.models](#)
[numpy](#)

[pylab](#)
[sys](#)

[unittest](#)

Functions

```
features(Signal=None, SamplingRate=1000.0, Filter={})
    Retrieves relevant EDA signal features.

    Kwargs:
        Signal (array): input EDA signal

        SamplingRate (float): sampling frequency (Hz)

        Filter (dict): Filter parameters

    Kwrvals:
        Signal (array): output filtered signal (see notes 1)

        Amplitude (array): signal pulses amplitudes (in the units of the input signal)

        Onset (array): indexes (or instants in seconds, see notes 2.a) of the SCRs onsets

        Peak (array): indexes (or instants in seconds, see notes 2.a) of the SCRs peaks

        Rise (array): SCRs rise times (in seconds)

        HalfRecovery (array): SCRs half-recovery times (in seconds)

        mean (float): mean

        std (float): standard deviation

        var (float): variance

        skew (ndarray): skewness

        kurtosis (array): kurtosis

        ad (float): absolute deviation

    Configurable fields:
    {"name": "eda.features", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter"], "outputs": ["Signal", "Amplitude", "Onset", "Peak", "Rise", "HalfRecovery", "mean", "std", "var", "skew", "kurtosis", "ad"]}

    See Also:
        filt
        models.basicSCR
        tls.statsf

    Notes:
        1 - If a filter is given as a parameter, then the returned keyworded values dict has a 'Signal' key.

        2 - If the sampling rate is defined, then:
            a) keys 'onset', and 'peak' are converted to instants of occurrence in seconds.

    Example:

    References:
        .. [1]
```

```
filt(Signal=None, SamplingRate=1000.0, UpperCutoff=0.25, Order=2)
    Filters an input EDA signal.

    If only input signal is provide, it returns the filtered EDA signal
    assuming a 1000Hz sampling frequency and a default low-pass filter
    with a cutoff frequency of 0.25Hz.

    Kwargs:
        Signal (array): input signal.

        SamplingRate (float): Sampling frequency (Hz).

        UpperCutoff (float): Low-pass filter cutoff frequency (Hz).

        LowerCutoff (float): High-pass filter cutoff frequency (Hz).

        Order (int): Filter order.

    Kwrvals:
        Signal (array): output filtered signal.
```

Configurable fields:{"name": "eda.filt", "config": {"UpperCutoff": "0.25", "SamplingRate": "1000.", "Order": "2"}, "inputs": ["Signal"], "out

See Also:
filt.zpdf

Notes:

Example:
Signal = load(...)
SamplingRate = ...
res = [filt](#)(Signal=Signal, SamplingRate=SamplingRate)
plot(res['Signal'])

References:
.. [1]

scl(Signal=None, SamplingRate=1000.0, Filter={'UpperCutoff': 0.05})
Kwargs:

Kwrvals:

Configurable fields:{"name": "eda.scl", "config": {"Filter": {"UpperCutoff": 0.05}, "SamplingRate": "1000.", "inputs": ["Signal"], "output:

See Also:

Notes:

Example:

References:
.. [1]

scr(Signal=None, SamplingRate=1000.0, Method='basic', Filter={})
Detects and extracts Skin Conductivity Responses (SCRs) information such as:
SCRs amplitudes, onsets, peak instant, rise, and half-recovery times.

Kwargs:
Signal (array): input EDA signal.
SamplingRate (float): Sampling frequency (Hz).
Method (string): SCR detection algorithm.
Filter (dict): filter parameters.

Kwrvals:
Signal (array): output filtered signal (see notes 1)
Amplitude (array): signal pulses amplitudes (in the units of the input signal)
Onset (array): indexes (or instants in seconds, see notes 2.a) of the SCRs onsets
Peak (array): indexes (or instants in seconds, see notes 2.a) of the SCRs peaks
Rise (array): SCRs rise times (in seconds)
HalfRecovery (array): SCRs half-recovery times (in seconds)

Configurable fields:{"name": "eda.scr", "config": {"SamplingRate": "1000.", "Method": "basic"}, "inputs": ["Signal", "Filter"], "outputs":

See Also:
filt
models.basicSCR
models.KBKSCR

Notes:
1 - If a filter is given as a parameter, then the returned keyworded values dict has a 'Signal' key.
2 - If the sampling rate is defined, then:
a) keys 'onset', and 'peak' are converted to instants of occurrence in seconds.

Example:

References:
.. [1]

Modules

[numpy](#)[pylab](#)[sys](#)[unittest](#)

Functions

KBKSCR(Signal=None, SamplingRate=1000.0)

Detects and extracts Skin Conductivity Responses (SCRs) information such as: SCRs amplitudes, onsets, peak instant, rise, and half-recovery times.

Kwargs:

Signal (array): input EDA signal.

SamplingRate (float): Sampling frequency (Hz).

Kwrvals:

Signal (array): output filtered signal (see notes 1)

Amplitude (array): signal pulses amplitudes (in the units of the input signal)

Onset (array): indexes (or instants in seconds, see notes 2.a) of the SCRs onsets

Peak (array): indexes (or instants in seconds, see notes 2.a) of the SCRs peaks

TODO: Rise (array): SCRs rise times (in seconds)

TODO: HalfRecovery (array): SCRs half-recovery times (in seconds)

See Also:

flt.zpdf

Notes:

1 - If the sampling rate is defined, then:

a) keys 'onset', and 'peak' are converted to instants of occurrence in seconds.

2- Less sensitive than Gamboa algorithm, but does not solve the overlapping SCRs problem.

Example:

References:

.. [1] K.H. Kim, S.W. Bang, and S.R. Kim
"Emotion recognition system using short-term monitoring of physiological signals"
Med. Biol. Eng. Comput., 2004, 42, 419-427

basicSCR(Signal=None, SamplingRate=1000.0, Filter={})

Detects and extracts Skin Conductivity Responses (SCRs) information such as: SCRs amplitudes, onsets, peak instant, rise, and half-recovery times.

Kwargs:

Signal (array): input EDA signal.

SamplingRate (float): Sampling frequency (Hz).

Method (string): SCR detection algorithm.

Filter (dict): filter parameters.

Kwrvals:

Signal (array): output filtered signal (see notes 1)

Amplitude (array): signal pulses amplitudes (in the units of the input signal)

Onset (array): indexes (or instants in seconds, see notes 2.a) of the SCRs onsets

Peak (array): indexes (or instants in seconds, see notes 2.a) of the SCRs peaks

Rise (array): SCRs rise times (in seconds)

HalfRecovery (array): SCRs half-recovery times (in seconds)

See Also:

`filt`

Notes:

1 - If a filter is given as a parameter, then the returned keyworded values dict has a 'Signal' key.

2 - If the sampling rate is defined, then:

a) keys 'onset', and 'peak' are converted to instants of occurrence in seconds.

Example:

References:

.. [1]

gamboa(to, tp, yo, yp)

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

EEG functions

Modules

[scipy.interpolate](#)
[numpy](#)[pandas](#)
[pylab](#)[scipy.signal](#)
[scipy.stats.stats](#)

Functions

EMDiit(Signal=None, niter=20, coef=0.6)**HSPlot**(T, W, matrix, NTicks=8)**HilbertSpectrum**(IMF, Fs, NF=100)**ICAFilter**(signal=None)**analyticSignal**(signal=None, axis=-1)**bandPower**(signal=None, Fs=1, lower=0, upper=1, NFFT=4096, axis=-1)**car**(signal=None)**emd**(data, extrapolation='mirror', nimfs=12, sifting_distance=0.2)

Perform a Empirical Mode Decomposition on a data set.

This function will return an array of all the Emperical Mode Functions as defined in [1]_, which can be used for further Hilbert Spectral Analysis.

The EMD uses a spline interpolation function to approximate the upper and lower envelopes of the signal, this routine implements a extrapolation routine as described in [2]_ as well as the standard spline routine. The extrapolation method removes the artifacts introduced by the spline fit at the ends of the data set, by making the dataset a continuous circle.

Reproduced from github.com/jaidevd/pyhht

Parameters

data : array_like
Signal Data

extrapolation : str, optional
Sets the extrapolation method for edge effects.
Options: None
 'mirror'
Default: 'mirror'

nimfs : int, optional
Sets the maximum number of IMFs to be found
Default : 12

sifting_distance : float, optional
Sets the minimum variance between IMF iterations.
Default : 0.2

Returns

IMFs : ndarray
An array of shape (len(data),N) where N is the number of found IMFs

Notes

References

.. [1] Huang H. et al. 1998 'The empirical mode decomposition and the Hilbert spectrum for nonlinear and non-stationary time series analysis.'
Proceedings of the Royal Society 454, 903-995

.. [2] Zhao J., Huang D. 2001 'Mirror extending and circular spline function for empirical mode decomposition method'
Journal of Zhejiang University (Science) V.2, No.3,P247-252

.. [3] Rato R.T., Ortigueira M.D., Batista A.G 2008 'On the HHT, its problems, and some solutions.'
Mechanical Systems and Signal Processing 22 1374-1394

intervalThresholding(IMF=None, thr=None)

laplacian(signal=None, indList=None)

noiseMeter(data=None, Fs=None, norm=False)

plf(signal=None, pairList=None)

power(signal=None, NFFT=4096, axis=-1)

spectPlot(X, Y, Z, **kwarg)

syncLikelihood(signal=None)

windower(signal=None, Fs=None, length=128, shift=64, fcn=None, fcnArgs={})


```
.. module:: emg
    :platform: Unix, Windows
    :synopsis: This module provides various functions to handle EMG signals.

.. moduleauthor:: Filipe Canento
```

Modules

[numpy](#) [scipy](#) [unittest](#)
[pylab](#) [sys](#)

Functions

emg(Signal=None, SamplingRate=1000.0, Filter=())
 EMG signal processing and feature extraction.

Kwargs:

- Signal (array): input signal.
- SamplingRate (float): Sampling frequency (Hz).
- Filter (dict): filter parameters.

Kwrvals:

- Signal (array): output filtered signal (see notes 1).
- onoff (array): indexes of EMG onsets

Configurable fields:{"name": "emg.emg", "config": {"SamplingRate": "1000.0"}, "inputs": ["Signal", "Filter"], "outputs": ["Signal", "onoff"]}

See Also:

- filt
- onoff

Notes:

Example:

References:

- .. [1]

features(Signal=None, SamplingRate=1000.0, Filter={})
 Retrieves relevant EMG signal features.

Kwargs:

- Signal (array): input signal.
- SamplingRate (float): Sampling frequency (Hz).
- Filter (dict): filter parameters.

Kwrvals:

- Signal (array): output filtered signal (see notes 1).
- onoff (array): indexes of EMG onsets
- mean (float): mean
- std (float): standard deviation
- var (float): variance
- skew (ndarray): skewness
- kurtosis (array): kurtosis
- ad (float): absolute deviation

Configurable fields:

{"name": "emg.features", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter"], "outputs": ["Signal", "onoff", "mean", "std", 'ad']}

See Also:

- filt
- onoff
- tls.statsf

Notes:

Example:

References:

- .. [1]

filt(Signal=None, SamplingRate=1000.0, UpperCutoff=None, LowerCutoff=100.0, Order=4.0)
Filters an input EMG signal.

If only input signal is provide, it returns the filtered EMG signal assuming a 1000Hz sampling frequency and a default high-pass filter with a cutoff frequency of 100Hz.

Kwargs:

Signal (array): input signal.

SamplingRate (float): Sampling frequency (Hz).

UpperCutoff (float): Low-pass filter cutoff frequency (Hz).

LowerCutoff (float): High-pass filter cutoff frequency (Hz).

Order (int): Filter order.

Kwrvals:

Signal (array): output filtered signal.

Configurable fields:

{"name": "emg.filt", "config": {"SamplingRate": "1000.", "LowerCutoff": "100.", "Order": "4."}, "inputs": ["Signal", "UpperCutoff"], "outputs": ["Signal"]}

See Also:

flt.zpdfr

Notes:

Example:

```
Signal = load(...)
SamplingRate = ...
res = filt(Signal=Signal, SamplingRate=SamplingRate)
plot(res['Signal'])
```

References:

.. [1]

onoff(Signal=None, SamplingRate=1000.0, Thres=None, ws=50.0)
EMG signal onset detection.

Kwargs:

Signal (array): input signal.

SamplingRate (float): Sampling frequency (Hz).

Thres (float): detection threshold

ws (float): detection window size in milliseconds

Kwrvals:

onoff (array): indexes of EMG onsets

Configurable fields:{"name": "emg.onoff", "config": {"SamplingRate": "1000.", "ws": "50.0"}, "inputs": ["Signal", "Thres"], "outputs": ["onoff"]}

See Also:

flt.smooth

Notes:

Example:

References:

.. [1]

Modules

[numpy](#)
[pylab](#)[scipy.signal](#)
[unittest](#)[scipy.signal.windows](#)

Functions

filterSignal(Signal=None, SamplingRate=1000.0, FilterType='FIR', Order=None, Frequency=None, BandType='lowpass', **kwargs)

Filter a signal according to the given parameters.

Uses a forward-backward filter implementation. Therefore, the combined filter has linear phase.

Supported filter functions: FIR, Butterworth, Chebyshev Type I, Chebyshev Type II, Elliptic, and Bessel.

Kwargs:

Signal (array): The signal to filter.

SamplingRate (int, float): The sampling frequency (Hz).

FilterType (str): The filter function: 'FIR', 'butter', 'cheby1', 'cheby2', 'ellip', or 'bessel' (default='FIR').

Order (int): Order of the filter.

Frequency (int, float, list, array): The cutoff frequency (or list/array of low and high cutoff frequencies).

BandType (str): The type of the filter: 'lowpass', 'highpass', 'bandpass', or 'bandstop' (default='lowpass').

**kwargs (dict): Additional keyword arguments are passed to the underlying scipy.signal function.

Kwrvals:

Signal (array): The filtered signal.

SamplingRate (float): The sampling frequency (Hz).

Filter (dict): The filter parameters.

See Also:

`_getFilter`

`plotFilter`

`scipy.signal.filtfilt`

Notes:

Example:

References:

.. [1]

firfilt(data, n, l, h, SamplingRate)

--DEPRECATED--

plotFilter(FilterType='FIR', Order=None, Frequency=None, SamplingRate=1000.0, BandType='lowpass', path=None, show=True, **kwargs)

Plot the frequency response of the filter specified with the given parameters.

Supported filter functions: FIR, Butterworth, Chebyshev Type I, Chebyshev Type II, Elliptic, and Bessel.

Kwargs:

FilterType (str): The filter function: 'FIR', 'butter', 'cheby1', 'cheby2', 'ellip', or 'bessel' (default='FIR').

Order (int): Order of the filter.

Frequency (int, float, list, array): The cutoff frequency (or array of frequencies).

SamplingRate (int, float): The sampling frequency (Hz).

BandType (str): The type of the filter: 'lowpass', 'highpass', 'bandpass', or 'bandstop' (default='lowpass').

path (str): If given, the plot will be saved to the file specified (default=None).

show (bool): If True, show the plot immediately.

****kwargs** (dict): Additional keyword arguments are passed to the underlying `scipy.signal` function.

Kwrvals:

See Also:

`filterSignal`

Notes:

Example:

References:

.. [1]

smooth(Signal=None, Window={})

Smooth data using a N-point moving average filter.

This implementation uses the convolution of a filter kernel with the input signal `x` to compute the smoothed signal.

Parameters

`x` : ndarray

Input signal data.

`n` : int, float, ndarray, optional

Number of points of the filter kernel. If this is an `int`, the filter kernel will have `n` points. If this is a `float`, the number of points for the filter kernel will be set as `n*size(x)`. If this is a `ndarray`, it will be directly taken as the filter kernel.

Default: 10

`wtype` : str, function, optional

Method that should be used to determine the filter kernel. If this is a `function`, it will be invoked with parameters `n` and `args` to determine the filter kernel. If this is a `str`, the function with the matching or most similar name belonging to the module `scipy.signal.windows` is used.

Default: `boxzen`

`*args` : optional

Additional parameters that may be required by the filter kernel function

Returns

`y` : ndarray

A smoothed version of the input signal computed using a filter kernel of size `n` generated according to `wtype`

See Also

`scipy.signal.windows`

Notes

A combination smoothing method `boxzen` was introduced and is currently used as default to produce the output signal.

This method first smooths the signal using the `scipy.windows.boxcar` window and then smooths it again using the `scipy.windows.parzen` window.

The resulting signals can be quite interesting, as `boxcar` retains a great proximity to the original data waveforms, and `parzen` removes the rough edges.

Example

```
t = arange(0,2*pi,.1)
x = sin(t)+0.5*(rand(len(t))-0.5)
y = smooth(x)
plot(x, '.')
plot(y)
legend(('original data', 'smoothed with boxzen'))
```

References

.. [1] Wikipedia, "Moving Average". http://en.wikipedia.org/wiki/Moving_average

.. [2] S. W. Smith, "Moving Average Filters - Implementation by Convolution". <http://www.dspguide.com/ch15/1.htm>

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

windowfcn(w)

Retrieve the appropriate window function that corresponds to the descriptor `w`.

Parameters

w : str, function
Name or descriptor of the window function. If this is a `function`, it will be directly returned. If this is a `str`, a window function with a matching name will be searched for in module ``scipy.signal.windows``.

Returns

f : function
The window function corresponding to the descriptor `w`

See Also

scipy.signal.windows

Notes

If no window function with a name matching the descriptor `w` is found in the module ``scipy.signal.windows``, the window function with the most similar name is be returned and a warning is issued.

Example

```
f = windowfcn('gaussian')  
f = windowfcn('gauss')
```

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

zpdfr(Signal=None, SamplingRate=None, UpperCutoff=None, LowerCutoff=None, Order=4.0)

--DEPRECATED FUNCTION--

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

Modules

[numpy](#)

Functions

hr(Signal=None, SamplingRate=1000.0)

ToDo.

Parameters

ToDo : ToDo

ToDo.

Returns

kwrvls : dict

A keyworded return values dict is returned with the following keys:

IBI : ndarray

ToDo.

HR : ndarray

ToDo.

See Also

bvp.pulse

ecg.ecg

Example

ToDo

References

.. ToDo

```
.. module:: peakd
  :platform: Unix, Windows
  :synopsis: This module provides various methods for peak detection.

.. moduleauthor:: Filipe Canento, Carlos Carreiras, Francisco David
```

Modules

[bspyclone.filt](#)[numpy](#)

Functions

sgndiff(Signal=None, a=-1)

Determines Signal peaks.

Kwargs:

Signal (array): input signal.

a (int): Whether to return maxima (a = -1, the default) or minima (a = 1) points.

Kwrvals:

Peak (array): peak indexes.

a is -1 as default, to detect peaks. if you want to detect minimum values, then change to 1

See Also:

Notes:

Example:

References:

.. [1]

ssf(Signal=None, SamplingRate=1000.0, Filter={})

Determines Signal peaks.

Kwargs:

Signal (array): input signal

SamplingRate (float): Sampling frequency (Hz)

Filter (dict): Filter coefficients

Kwrvals:

Signal (array):

Onset (array):

SSF (array):

See Also:

Notes:

Example:

References:

.. [1] W.Zong, T.Heldt, G.B. Moody, and R.G. Mark,
"An Open-
source Algorithm to Detect Onset of Arterial Blood Pressure Pulses",
Computers in Cardiology 2003; 30:259-262

This module provides various functions to ...

Functions:

[loadbpf\(\)](#)

Modules

[copy](#)

[numpy](#)

[pylab](#)

Classes

[numpy.ndarray](#)([__builtin__.object](#))

[bpararray](#)

class **bpararray**([numpy.ndarray](#))

#-----

Method resolution order:

[bpararray](#)

[numpy.ndarray](#)

[__builtin__.object](#)

Methods defined here:

__array_finalize__(self, obj)

toADC(self)

toV(self)

tomV(self)

touS(self)

Static methods defined here:

__new__(cls, ndarr, hdr={})

Data descriptors defined here:

__dict__
dictionary for instance variables (if defined)

Methods inherited from [numpy.ndarray](#):

__abs__(...)
x.[__abs__](#)() <=> abs(x)

__add__(...)
x.[__add__](#)(y) <=> x+y

__and__(...)

x.[__and__](#)(y) <==> x&y

[__array__](#)(...)

a.[__array__](#)(|dtype) -> reference if type unchanged, copy otherwise.

Returns either a new reference to self if dtype is not given or a new array of provided data type if dtype is different from the current dtype of the array.

[__array_prepare__](#)(...)

a.[__array_prepare__](#)(obj) -> Object of same type as [ndarray](#) object obj.

[__array_wrap__](#)(...)

a.[__array_wrap__](#)(obj) -> Object of same type as [ndarray](#) object a.

[__contains__](#)(...)

x.[__contains__](#)(y) <==> y in x

[__copy__](#)(...)

a.[__copy__](#)([order])

Return a copy of the array.

Parameters

order : {'C', 'F', 'A'}, optional

If order is 'C' (False) then the result is contiguous (default).

If order is 'Fortran' (True) then the result has fortran order.

If order is 'Any' (None) then the result has fortran order

only if the array already is in fortran order.

[__deepcopy__](#)(...)

a.[__deepcopy__](#)() -> Deep copy of array.

Used if copy.deepcopy is called on an array.

[__delitem__](#)(...)

x.[__delitem__](#)(y) <==> del x[y]

[__delslice__](#)(...)

x.[__delslice__](#)(i, j) <==> del x[i:j]

Use of negative indices is not supported.

[__div__](#)(...)

x.[__div__](#)(y) <==> x/y

[__divmod__](#)(...)

x.[__divmod__](#)(y) <==> divmod(x, y)

[__eq__](#)(...)

x.[__eq__](#)(y) <==> x==y

[__float__](#)(...)

x.[__float__](#)() <==> float(x)

[__floordiv__](#)(...)

x.[__floordiv__](#)(y) <==> x//y

[__ge__](#)(...)

x.[__ge__](#)(y) <==> x>=y

[__getitem__](#)(...)

x.[__getitem__](#)(y) <==> x[y]

[__getslice__](#)(...)

x.[__getslice__](#)(i, j) <==> x[i:j]

Use of negative indices is not supported.

__gt__(...)
x.**__gt__**(y) <==> x>y

__hex__(...)
x.**__hex__**() <==> hex(x)

__iadd__(...)
x.**__iadd__**(y) <==> x+=y

__iand__(...)
x.**__iand__**(y) <==> x&y

__idiv__(...)
x.**__idiv__**(y) <==> x/=y

__ifloordiv__(...)
x.**__ifloordiv__**(y) <==> x//y

__ilshift__(...)
x.**__ilshift__**(y) <==> x<<y

__imod__(...)
x.**__imod__**(y) <==> x%=y

__imul__(...)
x.**__imul__**(y) <==> x*=y

__index__(...)
x[y:z] <==> x[y.**__index__**():z.**__index__**()]

__int__(...)
x.**__int__**() <==> int(x)

__invert__(...)
x.**__invert__**() <==> ~x

__ior__(...)
x.**__ior__**(y) <==> x|=y

__ipow__(...)
x.**__ipow__**(y) <==> x**=y

__irshift__(...)
x.**__irshift__**(y) <==> x>>y

__isub__(...)
x.**__isub__**(y) <==> x-=y

__iter__(...)
x.**__iter__**() <==> iter(x)

__itruediv__(...)
x.**__itruediv__**(y) <==> x/y

__ixor__(...)
x.**__ixor__**(y) <==> x^=y

__le__(...)
x.**__le__**(y) <==> x<=y

__len__(...)
x.**__len__**() <==> len(x)

__long__(...)
x.**__long__**() <==> long(x)

__lshift__(...)
x.**__lshift__**(y) <==> x<<y

```

__lt__(...)
    x.__lt__(y) <==> x<y

__mod__(...)
    x.__mod__(y) <==> x%y

__mul__(...)
    x.__mul__(y) <==> x*y

__ne__(...)
    x.__ne__(y) <==> x!=y

__neg__(...)
    x.__neg__() <==> -x

__nonzero__(...)
    x.__nonzero__() <==> x != 0

__oct__(...)
    x.__oct__() <==> oct(x)

__or__(...)
    x.__or__(y) <==> x|y

__pos__(...)
    x.__pos__() <==> +x

__pow__(...)
    x.__pow__(y[, z]) <==> pow(x, y[, z])

__radd__(...)
    x.__radd__(y) <==> y+x

__rand__(...)
    x.__rand__(y) <==> y&x

__rdiv__(...)
    x.__rdiv__(y) <==> y/x

__rdivmod__(...)
    x.__rdivmod__(y) <==> divmod(y, x)

__reduce__(...)
    a.__reduce__()

    For pickling.

__repr__(...)
    x.__repr__() <==> repr(x)

__rfloordiv__(...)
    x.__rfloordiv__(y) <==> y//x

__rlshift__(...)
    x.__rlshift__(y) <==> y<<x

__rmod__(...)
    x.__rmod__(y) <==> y%x

__rmul__(...)
    x.__rmul__(y) <==> y*x

__ror__(...)
    x.__ror__(y) <==> y|x

__rpow__(...)
    y.__rpow__(x[, z]) <==> pow(x, y[, z])

```

```

__rrshift__(...)
    x.\_\_rrshift\_\_(y) <==> y>>x

__rshift__(...)
    x.\_\_rshift\_\_(y) <==> x>>y

__rsub__(...)
    x.\_\_rsub\_\_(y) <==> y-x

__rtruediv__(...)
    x.\_\_rtruediv\_\_(y) <==> y/x

__rxor__(...)
    x.\_\_rxor\_\_(y) <==> y^x

__setitem__(...)
    x.\_\_setitem\_\_(i, y) <==> x[i]=y

__setslice__(...)
    x.\_\_setslice\_\_(i, j, y) <==> x[i:j]=y

    Use of negative indices is not supported.

__setstate__(...)
    a.\_\_setstate\_\_(version, shape, dtype, isfortran, rawdata)

    For unpickling.

    Parameters
    -----
    version : int
        optional pickle version. If omitted defaults to 0.
    shape : tuple
    dtype : data-type
    isFortran : bool
    rawdata : string or list
        a binary string with the data (or a list if 'a' is an object array)

__str__(...)
    x.\_\_str\_\_() <==> str(x)

__sub__(...)
    x.\_\_sub\_\_(y) <==> x-y

__truediv__(...)
    x.\_\_truediv\_\_(y) <==> x/y

__xor__(...)
    x.\_\_xor\_\_(y) <==> x^y

all(...)
    a.all(axis=None, out=None)

    Returns True if all elements evaluate to True.

    Refer to `numpy.all` for full documentation.

    See Also
    -----
    numpy.all : equivalent function

any(...)
    a.any(axis=None, out=None)

    Returns True if any of the elements of `a` evaluate to True.

    Refer to `numpy.any` for full documentation.

    See Also
    -----
    numpy.any : equivalent function

```

argmax(...)

a.[argmax](#)(axis=None, out=None)

Return indices of the maximum values along the given axis.

Refer to ``numpy.argmax`` for full documentation.

See Also

`numpy.argmax` : equivalent function

argmin(...)

a.[argmin](#)(axis=None, out=None)

Return indices of the minimum values along the given axis of ``a``.

Refer to ``numpy.argmin`` for detailed documentation.

See Also

`numpy.argmin` : equivalent function

argpartition(...)

a.[argpartition](#)(kth, axis=-1, kind='introselect', order=None)

Returns the indices that would partition this array.

Refer to ``numpy.argpartition`` for full documentation.

.. versionadded:: 1.8.0

See Also

`numpy.argpartition` : equivalent function

argsort(...)

a.[argsort](#)(axis=-1, kind='quicksort', order=None)

Returns the indices that would sort this array.

Refer to ``numpy.argsort`` for full documentation.

See Also

`numpy.argsort` : equivalent function

astype(...)

a.[astype](#)(dtype, order='K', casting='unsafe', subok=True, copy=True)

Copy of the array, cast to a specified type.

Parameters

`dtype` : str or dtype

Typecode or data-type to which the array is cast.

`order` : {'C', 'F', 'A', 'K'}, optional

Controls the memory layout order of the result.

'C' means C order, 'F' means Fortran order, 'A'

means 'F' order if all the arrays are Fortran contiguous,

'C' order otherwise, and 'K' means as close to the

order the array elements appear in memory as possible.

Default is 'K'.

`casting` : {'no', 'equiv', 'safe', 'same_kind', 'unsafe'}, optional

Controls what kind of data casting may occur. Defaults to 'unsafe' for backwards compatibility.

* 'no' means the data types should not be cast at all.

* 'equiv' means only byte-order changes are allowed.

* 'safe' means only casts which can preserve values are allowed.

* 'same_kind' means only safe casts or casts within a kind, like float64 to float32, are allowed.

* 'unsafe' means any data conversions may be done.

subok : bool, optional
If True, then sub-classes will be passed-through (default), otherwise the returned array will be forced to be a base-class array.
copy : bool, optional
By default, `astype` always returns a newly allocated array. If this is set to false, and the ``dtype``, ``order``, and ``subok`` requirements are satisfied, the input array is returned instead of a copy.

Returns

`arr_t` : [ndarray](#)
Unless ``copy`` is False and the other conditions for returning the input array are satisfied (see description for ``copy`` input parameter), ``arr_t`` is a new array of the same shape as the input array, with dtype, order given by ``dtype``, ``order``.

Notes

Starting in NumPy 1.9, `astype` method now returns an error if the string dtype to cast to is not long enough in 'safe' casting mode to hold the max value of integer/float array that is being casted. Previously the casting was allowed even if the result was truncated.

Raises

ComplexWarning
When casting from complex to float or int. To avoid this, one should use ``a.real.astype(t)``.

Examples

>>> x = np.array([1, 2, 2.5])
>>> x
array([1. , 2. , 2.5])

>>> x.[astype](#)(int)
array([1, 2, 2])

byteswap(...)

`a.byteswap(inplace)`

Swap the bytes of the array elements

Toggle between low-endian and big-endian data representation by returning a byteswapped array, optionally swapped in-place.

Parameters

`inplace` : bool, optional
If ``True``, swap bytes in-place, default is ``False``.

Returns

`out` : [ndarray](#)
The byteswapped array. If ``inplace`` is ``True``, this is a view to self.

Examples

>>> A = np.array([1, 256, 8755], dtype=np.int16)
>>> map(hex, A)
['0x1', '0x100', '0x2233']
>>> A.[byteswap](#)(True)
array([256, 1, 13090], dtype=int16)
>>> map(hex, A)
['0x100', '0x1', '0x3322']

Arrays of strings are not swapped

>>> A = np.array(['ceg', 'fac'])
>>> A.[byteswap](#)()
array(['ceg', 'fac'],
 dtype='<S3')

choose(...)

a.[choose](#)(choices, out=None, mode='raise')

Use an index array to construct a new array from a set of choices.

Refer to ``numpy.choose`` for full documentation.

See Also

`numpy.choose` : equivalent function

clip(...)

a.[clip](#)(a_min, a_max, out=None)

Return an array whose values are limited to ```[a_min, a_max]```.

Refer to ``numpy.clip`` for full documentation.

See Also

`numpy.clip` : equivalent function

compress(...)

a.[compress](#)(condition, axis=None, out=None)

Return selected slices of this array along given axis.

Refer to ``numpy.compress`` for full documentation.

See Also

`numpy.compress` : equivalent function

conj(...)

a.[conj](#)()

Complex-conjugate all elements.

Refer to ``numpy.conjugate`` for full documentation.

See Also

`numpy.conjugate` : equivalent function

conjugate(...)

a.[conjugate](#)()

Return the complex conjugate, element-wise.

Refer to ``numpy.conjugate`` for full documentation.

See Also

`numpy.conjugate` : equivalent function

copy(...)

a.[copy](#)(order='C')

Return a copy of the array.

Parameters

order : {'C', 'F', 'A', 'K'}, optional

Controls the memory layout of the copy. 'C' means C-order, 'F' means F-order, 'A' means 'F' if ``a`` is Fortran contiguous, 'C' otherwise. 'K' means match the layout of ``a`` as closely as possible. (Note that this function and `:func:numpy.copy` are very similar, but have different default values for their order= arguments.)

See also

`numpy.copy`

`numpy.copyto`

Examples

```
>>> x = np.array([[1,2,3],[4,5,6]], order='F')
```

```
>>> y = x.copy()
```

```
>>> x.fill(0)
```

```
>>> x
array([[0, 0, 0],
       [0, 0, 0]])
```

```
>>> y
array([[1, 2, 3],
       [4, 5, 6]])
```

```
>>> y.flags['C_CONTIGUOUS']
True
```

cumprod(...)

a.[cumprod](#)(axis=None, dtype=None, out=None)

Return the cumulative product of the elements along the given axis.

Refer to ``numpy.cumprod`` for full documentation.

See Also

`numpy.cumprod` : equivalent function

cumsum(...)

a.[cumsum](#)(axis=None, dtype=None, out=None)

Return the cumulative sum of the elements along the given axis.

Refer to ``numpy.cumsum`` for full documentation.

See Also

`numpy.cumsum` : equivalent function

diagonal(...)

a.[diagonal](#)(offset=0, axis1=0, axis2=1)

Return specified diagonals. In NumPy 1.9 the returned array is a read-only view instead of a copy as in previous NumPy versions. In NumPy 1.10 the read-only restriction will be removed.

Refer to `:func:`numpy.diagonal`` for full documentation.

See Also

`numpy.diagonal` : equivalent function

dot(...)

a.[dot](#)(b, out=None)

Dot product of two arrays.

Refer to ``numpy.dot`` for full documentation.

See Also

`numpy.dot` : equivalent function

Examples

```
>>> a = np.eye(2)
>>> b = np.ones((2, 2)) * 2
>>> a.dot(b)
array([[ 2.,  2.]])
```

```
[ 2.,  2.]])
```

This array method can be conveniently chained:

```
>>> a.dot(b).dot(b)
array([[ 8.,  8.],
       [ 8.,  8.]])
```

dump(...)

[a.dump](#)(file)

Dump a pickle of the array to the specified file.
The array can be read back with `pickle.load` or `numpy.load`.

Parameters

file : str

A string naming the dump file.

dumps(...)

[a.dumps](#)()

Returns the pickle of the array as a string.
`pickle.loads` or `numpy.loads` will convert the string back to an array.

Parameters

None

fill(...)

[a.fill](#)(value)

Fill the array with a scalar value.

Parameters

value : scalar

All elements of `a` will be assigned this value.

Examples

```
>>> a = np.array([1, 2])
>>> a.fill(0)
>>> a
array([0, 0])
>>> a = np.empty(2)
>>> a.fill(1)
>>> a
array([ 1.,  1.]])
```

flatten(...)

[a.flatten](#)(order='C')

Return a copy of the array collapsed into one dimension.

Parameters

order : {'C', 'F', 'A'}, optional

Whether to flatten in C (row-major), Fortran (column-major) order,
or preserve the C/Fortran ordering from `a`.
The default is 'C'.

Returns

y : [ndarray](#)

A copy of the input array, flattened to one dimension.

See Also

ravel : Return a flattened array.

flat : A 1-D flat iterator over the array.

Examples

```

-----
>>> a = np.array([[1,2], [3,4]])
>>> a.flatten()
array([1, 2, 3, 4])
>>> a.flatten('F')
array([1, 3, 2, 4])

```

getfield(...)

a.[getfield](#)(dtype, offset=0)

Returns a field of the given array as a certain type.

A field is a view of the array data with a given data-type. The values in the view are determined by the given type and the offset into the current array in bytes. The offset needs to be such that the view dtype fits in the array dtype; for example an array of dtype complex128 has 16-byte elements. If taking a view with a 32-bit integer (4 bytes), the offset needs to be between 0 and 12 bytes.

Parameters

dtype : str or dtype

The data type of the view. The dtype size of the view can not be larger than that of the array itself.

offset : int

Number of bytes to skip before beginning the element view.

Examples

```

>>> x = np.diag([1.+1.j]*2)
>>> x[1, 1] = 2 + 4.j
>>> x
array([[ 1.+1.j,  0.+0.j],
       [ 0.+0.j,  2.+4.j]])
>>> x.getfield(np.float64)
array([[ 1.,  0.],
       [ 0.,  2.]])

```

By choosing an offset of 8 bytes we can select the complex part of the array for our view:

```

>>> x.getfield(np.float64, offset=8)
array([[ 1.,  0.],
       [ 0.,  4.]])

```

item(...)

a.[item](#)(*args)

Copy an element of an array to a standard Python scalar and return it.

Parameters

*args : Arguments (variable number and type)

* none: in this case, the method only works for arrays with one element (``a.size == 1``), which element is copied into a standard Python scalar object and returned.

* int_type: this argument is interpreted as a flat index into the array, specifying which element to copy and return.

* tuple of int_types: functions as does a single int_type argument, except that the argument is interpreted as an nd-index into the array.

Returns

z : Standard Python scalar object

A copy of the specified element of the array as a suitable Python scalar

Notes

When the data type of ``a`` is longdouble or clongdouble, [item](#)() returns

a scalar array object because there is no available Python scalar that would not lose information. Void arrays return a buffer object for [item\(\)](#), unless fields are defined, in which case a tuple is returned.

``item`` is very similar to `a[args]`, except, instead of an array scalar, a standard Python scalar is returned. This can be useful for speeding up access to elements of the array and doing arithmetic on elements of the array using Python's optimized math.

Examples

```
-----
>>> x = np.random.randint(9, size=(3, 3))
>>> x
array([[3, 1, 7],
       [2, 8, 3],
       [8, 5, 3]])
>>> x.item(3)
2
>>> x.item(7)
5
>>> x.item((0, 1))
1
>>> x.item((2, 2))
3
```

itemset(...)

`a.itemset(*args)`

Insert scalar into an array (scalar is cast to array's dtype, if possible)

There must be at least 1 argument, and define the last argument as `*item`. Then, ``a.itemset(*args)`` is equivalent to but faster than ``a[args] = item``. The item should be a scalar value and ``args`` must select a single item in the array ``a``.

Parameters

```
-----
\*args : Arguments
    If one argument: a scalar, only used in case `a` is of size 1.
    If two arguments: the last argument is the value to be set
    and must be a scalar, the first argument specifies a single array
    element location. It is either an int or a tuple.
```

Notes

```
-----
Compared to indexing syntax, `itemset` provides some speed increase for placing a scalar into a particular location in an `ndarray`, if you must do this. However, generally this is discouraged: among other problems, it complicates the appearance of the code. Also, when using `itemset` (and `item`) inside a loop, be sure to assign the methods to a local variable to avoid the attribute look-up at each loop iteration.
```

Examples

```
-----
>>> x = np.random.randint(9, size=(3, 3))
>>> x
array([[3, 1, 7],
       [2, 8, 3],
       [8, 5, 3]])
>>> x.itemset(4, 0)
>>> x.itemset((2, 2), 9)
>>> x
array([[3, 1, 7],
       [2, 0, 3],
       [8, 5, 9]])
```

max(...)

`a.max(axis=None, out=None)`

Return the maximum along a given axis.

Refer to ``numpy.amax`` for full documentation.

See Also

numpy.amax : equivalent function

mean(...)

a.[mean](#)(axis=None, dtype=None, out=None)

Returns the average of the array elements along given axis.

Refer to ``numpy.mean`` for full documentation.

See Also

numpy.mean : equivalent function

min(...)

a.[min](#)(axis=None, out=None)

Return the minimum along a given axis.

Refer to ``numpy.amin`` for full documentation.

See Also

numpy.amin : equivalent function

newbyteorder(...)

arr.[newbyteorder](#)(new_order='S')

Return the array with the same data viewed with a different byte order.

Equivalent to::

```
arr.view(arr.dtype.newbytorder(new_order))
```

Changes are also made in all fields and sub-arrays of the array data type.

Parameters

new_order : string, optional

Byte order to force; a value from the byte order specifications above. ``new_order`` codes can be any of::

- * 'S' - swap dtype from current to opposite endian
- * {'<', 'L'} - little endian
- * {'>', 'B'} - big endian
- * {'=', 'N'} - native order
- * {'|', 'I'} - ignore (no change to byte order)

The default value ('S') results in swapping the current byte order. The code does a case-insensitive check on the first letter of ``new_order`` for the alternatives above. For example, any of 'B' or 'b' or 'bigish' are valid to specify big-endian.

Returns

new_arr : array

New array object with the dtype reflecting given change to the byte order.

nonzero(...)

a.[nonzero](#)()

Return the indices of the elements that are non-zero.

Refer to ``numpy.nonzero`` for full documentation.

See Also

numpy.nonzero : equivalent function

partition(...)

a.[partition](#)(kth, axis=-1, kind='introselect', order=None)

Rearranges the elements in the array in such a way that value of the element in kth position is in the position it would be in a sorted array. All elements smaller than the kth element are moved before this element and all equal or greater are moved behind it. The ordering of the elements in the two partitions is undefined.

.. versionadded:: 1.8.0

Parameters

kth : int or sequence of ints

Element index to partition by. The kth element value will be in its final sorted position and all smaller elements will be moved before it and all equal or greater elements behind it.

The order all elements in the partitions is undefined.

If provided with a sequence of kth it will partition all elements indexed by kth of them into their sorted position at once.

axis : int, optional

Axis along which to sort. Default is -1, which means sort along the last axis.

kind : {'introselect'}, optional

Selection algorithm. Default is 'introselect'.

order : list, optional

When `a` is an array with fields defined, this argument specifies which fields to compare first, second, etc. Not all fields need be specified.

See Also

numpy.partition : Return a partitioned copy of an array.

argpartition : Indirect partition.

sort : Full sort.

Notes

See ``np.partition`` for notes on the different algorithms.

Examples

```
>>> a = np.array([3, 4, 2, 1])
>>> a.partition(a, 3)
>>> a
array([2, 1, 3, 4])
```

```
>>> a.partition((1, 3))
array([1, 2, 3, 4])
```

prod(...)

a.[prod](#)(axis=None, dtype=None, out=None)

Return the product of the array elements over the given axis

Refer to `numpy.prod` for full documentation.

See Also

numpy.prod : equivalent function

ptp(...)

a.[ptp](#)(axis=None, out=None)

Peak to peak (maximum - minimum) value along a given axis.

Refer to `numpy.ptp` for full documentation.

See Also

numpy.ptp : equivalent function

put(...)

a.[put](#)(indices, values, mode='raise')

Set ``a.flat[n] = values[n]`` for all ``n`` in indices.

Refer to ``numpy.put`` for full documentation.

See Also

numpy.put : equivalent function

ravel(...)

a.[ravel](#)([order])

Return a flattened array.

Refer to ``numpy.ravel`` for full documentation.

See Also

numpy.ravel : equivalent function

[ndarray.flat](#) : a flat iterator on the array.

repeat(...)

a.[repeat](#)(repeats, axis=None)

Repeat elements of an array.

Refer to ``numpy.repeat`` for full documentation.

See Also

numpy.repeat : equivalent function

reshape(...)

a.[reshape](#)(shape, order='C')

Returns an array containing the same data with a new shape.

Refer to ``numpy.reshape`` for full documentation.

See Also

numpy.reshape : equivalent function

resize(...)

a.[resize](#)(new_shape, refcheck=True)

Change shape and size of array in-place.

Parameters

new_shape : tuple of ints, or ``n`` ints

Shape of resized array.

refcheck : bool, optional

If False, reference count will not be checked. Default is True.

Returns

None

Raises

ValueError

If ``a`` does not own its own data or references or views to it exist,
and the data memory must be changed.

SystemError

If the ``order`` keyword argument is specified. This behaviour is a
bug in NumPy.

See Also

resize : Return a new array with the specified shape.

Notes

This reallocates space for the data area if necessary.

Only contiguous arrays (data elements consecutive in memory) can be resized.

The purpose of the reference count check is to make sure you do not use this array as a buffer for another Python object and then reallocate the memory. However, reference counts can increase in other ways so if you are sure that you have not shared the memory for this array with another Python object, then you may safely set ``refcheck`` to False.

Examples

Shrinking an array: array is flattened (in the order that the data are stored in memory), resized, and reshaped:

```
>>> a = np.array([[0, 1], [2, 3]], order='C')
>>> a.resize((2, 1))
>>> a
array([[0],
       [1]])
```

```
>>> a = np.array([[0, 1], [2, 3]], order='F')
>>> a.resize((2, 1))
>>> a
array([[0],
       [2]])
```

Enlarging an array: as above, but missing entries are filled with zeros:

```
>>> b = np.array([[0, 1], [2, 3]])
>>> b.resize(2, 3) # new_shape parameter doesn't have to be a tuple
>>> b
array([[0, 1, 2],
       [3, 0, 0]])
```

Referencing an array prevents resizing...

```
>>> c = a
>>> a.resize((1, 1))
Traceback (most recent call last):
...
ValueError: cannot resize an array that has been referenced ...
```

Unless ``refcheck`` is False:

```
>>> a.resize((1, 1), refcheck=False)
>>> a
array([[0]])
>>> c
array([[0]])
```

round(...)

`a.round(decimals=0, out=None)`

Return ``a`` with each element rounded to the given number of decimals.

Refer to ``numpy.around`` for full documentation.

See Also

`numpy.around` : equivalent function

searchsorted(...)

`a.searchsorted(v, side='left', sorter=None)`

Find indices where elements of `v` should be inserted in `a` to maintain order.

For full documentation, see ``numpy.searchsorted``

See Also

`numpy.searchsorted` : equivalent function

setfield(...)

a.[setfield](#)(val, dtype, offset=0)

Put a value into a specified place in a field defined by a data-type.

Place ``val`` into ``a``'s field defined by ``dtype`` and beginning ``offset`` bytes into the field.

Parameters

val : object

Value to be placed in field.

dtype : dtype object

Data-type of the field in which to place ``val``.

offset : int, optional

The number of bytes into the field at which to place ``val``.

Returns

None

See Also

[getfield](#)

Examples

```
>>> x = np.eye(3)
>>> x.getfield(np.float64)
array([[ 1.,  0.,  0.],
       [ 0.,  1.,  0.],
       [ 0.,  0.,  1.]])
>>> x.setfield(3, np.int32)
>>> x.getfield(np.int32)
array([[3, 3, 3],
       [3, 3, 3],
       [3, 3, 3]])
>>> x
array([[ 1.00000000e+000,  1.48219694e-323,  1.48219694e-323],
       [ 1.48219694e-323,  1.00000000e+000,  1.48219694e-323],
       [ 1.48219694e-323,  1.48219694e-323,  1.00000000e+000]])
>>> x.setfield(np.eye(3), np.int32)
>>> x
array([[ 1.,  0.,  0.],
       [ 0.,  1.,  0.],
       [ 0.,  0.,  1.]])
```

setflags(...)

a.[setflags](#)(write=None, align=None, uic=None)

Set array flags WRITEABLE, ALIGNED, and UPDATEIFCOPY, respectively.

These Boolean-valued flags affect how numpy interprets the memory area used by ``a`` (see Notes below). The ALIGNED flag can only be set to True if the data is actually aligned according to the type. The UPDATEIFCOPY flag can never be set to True. The flag WRITEABLE can only be set to True if the array owns its own memory, or the ultimate owner of the memory exposes a writeable buffer interface, or is a string. (The exception for string is made so that unpickling can be done without copying memory.)

Parameters

write : bool, optional

Describes whether or not ``a`` can be written to.

align : bool, optional

Describes whether or not ``a`` is aligned properly for its type.

uic : bool, optional

Describes whether or not ``a`` is a copy of another "base" array.

Notes

Array flags provide information about how the memory area used for the array is to be interpreted. There are 6 Boolean flags in use, only three of which can be changed by the user: `UPDATEIFCOPY`, `WRITEABLE`, and `ALIGNED`.

`WRITEABLE` (W) the data area can be written to;

`ALIGNED` (A) the data and strides are aligned appropriately for the hardware (as determined by the compiler);

`UPDATEIFCOPY` (U) this array is a copy of some other array (referenced by `.base`). When this array is deallocated, the base array will be updated with the contents of this array.

All flags can be accessed using their first (upper case) letter as well as the full name.

Examples

```
>>> y
array([[3, 1, 7],
       [2, 0, 0],
       [8, 5, 9]])
>>> y.flags
  C_CONTIGUOUS : True
  F_CONTIGUOUS : False
  OWNDATA : True
  WRITEABLE : True
  ALIGNED : True
  UPDATEIFCOPY : False
>>> y.setflags(write=0, align=0)
>>> y.flags
  C_CONTIGUOUS : True
  F_CONTIGUOUS : False
  OWNDATA : True
  WRITEABLE : False
  ALIGNED : False
  UPDATEIFCOPY : False
>>> y.setflags(uic=1)
Traceback (most recent call last):
  File "<stdin>", line 1, in <module>
ValueError: cannot set UPDATEIFCOPY flag to True
```

sort(...)

`a.sort(axis=-1, kind='quicksort', order=None)`

Sort an array, in-place.

Parameters

`axis` : int, optional
Axis along which to sort. Default is -1, which means sort along the last axis.

`kind` : {'quicksort', 'mergesort', 'heapsort'}, optional
Sorting algorithm. Default is 'quicksort'.

`order` : list, optional
When ``a`` is an array with fields defined, this argument specifies which fields to compare first, second, etc. Not all fields need be specified.

See Also

`numpy.sort` : Return a sorted copy of an array.
`argsort` : Indirect sort.
`lexsort` : Indirect stable sort on multiple keys.
`searchsorted` : Find elements in sorted array.
`partition` : Partial sort.

Notes

See ```sort``` for notes on the different sorting algorithms.

Examples

```
>>> a = np.array([[1,4], [3,1]])
>>> a.sort(axis=1)
>>> a
array([[1, 4],
       [1, 3]])
>>> a.sort(axis=0)
>>> a
array([[1, 3],
       [1, 4]])
```

Use the ``order`` keyword to specify a field to use when sorting a structured array:

```
>>> a = np.array([('a', 2), ('c', 1)], dtype=[('x', 'S1'), ('y', int)])
>>> a.sort(order='y')
>>> a
array([('c', 1), ('a', 2)],
      dtype=[('x', '|S1'), ('y', '<i4')])
```

squeeze(...)

`a.squeeze(axis=None)`

Remove single-dimensional entries from the shape of ``a``.

Refer to ``numpy.squeeze`` for full documentation.

See Also

`numpy.squeeze` : equivalent function

std(...)

`a.std(axis=None, dtype=None, out=None, ddof=0)`

Returns the standard deviation of the array elements along given axis.

Refer to ``numpy.std`` for full documentation.

See Also

`numpy.std` : equivalent function

sum(...)

`a.sum(axis=None, dtype=None, out=None)`

Return the sum of the array elements over the given axis.

Refer to ``numpy.sum`` for full documentation.

See Also

`numpy.sum` : equivalent function

swapaxes(...)

`a.swapaxes(axis1, axis2)`

Return a view of the array with ``axis1`` and ``axis2`` interchanged.

Refer to ``numpy.swapaxes`` for full documentation.

See Also

`numpy.swapaxes` : equivalent function

take(...)

`a.take(indices, axis=None, out=None, mode='raise')`

Return an array formed from the elements of ``a`` at the given indices.

Refer to ``numpy.take`` for full documentation.

See Also

numpy.take : equivalent function

tobytes(...)

a.[tobytes](#)(order='C')

Construct Python bytes containing the raw data bytes in the array.

Constructs Python bytes showing a copy of the raw contents of data memory. The bytes object can be produced in either 'C' or 'Fortran', or 'Any' order (the default is 'C'-order). 'Any' order means C-order unless the F_CONTIGUOUS flag in the array is set, in which case it means 'Fortran' order.

.. versionadded:: 1.9.0

Parameters

order : {'C', 'F', None}, optional
Order of the data for multidimensional arrays:
C, Fortran, or the same as for the original array.

Returns

s : bytes
Python bytes exhibiting a copy of `a`'s raw data.

Examples

>>> x = np.array([[0, 1], [2, 3]])
>>> x.[tobytes](#)()
b'\x00\x00\x00\x00\x01\x00\x00\x00\x02\x00\x00\x00\x03\x00\x00\x00'
>>> x.[tobytes](#)('C') == x.[tobytes](#)()
True
>>> x.[tobytes](#)('F')
b'\x00\x00\x00\x00\x02\x00\x00\x00\x01\x00\x00\x00\x03\x00\x00\x00'

tofile(...)

a.[tofile](#)(fid, sep="", format="%s")

Write array to a file as text or binary (default).

Data is always written in 'C' order, independent of the order of `a`. The data produced by this method can be recovered using the function `fromfile()`.

Parameters

fid : file or str
An open file object, or a string containing a filename.
sep : str
Separator between array items for text output.
If "" (empty), a binary file is written, equivalent to
``file.write(a.[tobytes](#)())``.
format : str
Format string for text file output.
Each entry in the array is formatted to text by first converting it to the closest Python type, and then using "format" % item.

Notes

This is a convenience function for quick storage of array data. Information on endianness and precision is lost, so this method is not a good choice for files intended to archive data or transport data between machines with different endianness. Some of these problems can be overcome by outputting the data as text files, at the expense of speed and file size.

tolist(...)

a.[tolist](#)()

Return the array as a (possibly nested) list.

Return a copy of the array data as a (nested) Python list.
Data items are converted to the nearest compatible Python type.

Parameters

none

Returns

y : list

The possibly nested list of array elements.

Notes

The array may be recreated, ``a = np.array(a.[tolist\(\)](#))``.

Examples

```
>>> a = np.array([1, 2])
>>> a.tolist\(\)
[1, 2]
>>> a = np.array([[1, 2], [3, 4]])
>>> list(a)
[array([1, 2]), array([3, 4])]
>>> a.tolist\(\)
[[1, 2], [3, 4]]
```

tostring(...)

a.[tostring](#)(order='C')

Construct Python bytes containing the raw data bytes in the array.

Constructs Python bytes showing a copy of the raw contents of data memory. The bytes object can be produced in either 'C' or 'Fortran', or 'Any' order (the default is 'C'-order). 'Any' order means C-order unless the F_CONTIGUOUS flag in the array is set, in which case it means 'Fortran' order.

This function is a compatibility alias for tobytes. Despite its name it returns bytes not strings.

Parameters

order : {'C', 'F', None}, optional

Order of the data for multidimensional arrays:

C, Fortran, or the same as for the original array.

Returns

s : bytes

Python bytes exhibiting a copy of `a`'s raw data.

Examples

```
>>> x = np.array([[0, 1], [2, 3]])
>>> x.tobytes\(\)
b'\x00\x00\x00\x00\x01\x00\x00\x00\x02\x00\x00\x00\x03\x00\x00\x00'
>>> x.tobytes('C') == x.tobytes()
True
>>> x.tobytes('F')
b'\x00\x00\x00\x00\x02\x00\x00\x00\x01\x00\x00\x00\x03\x00\x00\x00'
```

trace(...)

a.[trace](#)(offset=0, axis1=0, axis2=1, dtype=None, out=None)

Return the sum along diagonals of the array.

Refer to `numpy.trace` for full documentation.

See Also

numpy.trace : equivalent function

transpose(...)

a.[transpose](#)(*axes)

Returns a view of the array with axes transposed.

For a 1-D array, this has no effect. (To change between column and row vectors, first cast the 1-D array into a matrix object.)

For a 2-D array, this is the usual matrix transpose.

For an n-D array, if axes are given, their order indicates how the axes are permuted (see Examples). If axes are not provided and

``a.shape = (i[0], i[1], ... i[n-2], i[n-1])``, then
``a.transpose()``.shape = (i[n-1], i[n-2], ... i[1], i[0])`.

Parameters

axes : None, tuple of ints, or `n` ints

- * None or no argument: reverses the order of the axes.
- * tuple of ints: `i` in the `j`-th place in the tuple means `a`'s `i`-th axis becomes `a.[transpose](#)()`'s `j`-th axis.
- * `n` ints: same as an n-tuple of the same ints (this form is intended simply as a "convenience" alternative to the tuple form)

Returns

out : [ndarray](#)

View of `a`, with axes suitably permuted.

See Also

[ndarray.T](#) : Array property returning the array transposed.

Examples

```
>>> a = np.array([[1, 2], [3, 4]])
>>> a
array([[1, 2],
       [3, 4]])
>>> a.transpose()
array([[1, 3],
       [2, 4]])
>>> a.transpose((1, 0))
array([[1, 3],
       [2, 4]])
>>> a.transpose(1, 0)
array([[1, 3],
       [2, 4]])
```

var(...)

a.[var](#)(axis=None, dtype=None, out=None, ddof=0)

Returns the variance of the array elements, along given axis.

Refer to `numpy.var` for full documentation.

See Also

`numpy.var` : equivalent function

view(...)

a.[view](#)(dtype=None, type=None)

New view of array with the same data.

Parameters

dtype : data-type or [ndarray](#) sub-class, optional

Data-type descriptor of the returned view, e.g., float32 or int16. The default, None, results in the view having the same data-type as `a`.

This argument can also be specified as an [ndarray](#) sub-class, which then specifies the type of the returned object (this is equivalent to

setting the ``type`` parameter).
type : Python type, optional
Type of the returned view, e.g., [ndarray](#) or matrix. Again, the default None results in type preservation.

Notes

``a.view()`` is used two different ways:

``a.view(some_dtype)`` or ``a.view(dtype=some_dtype)`` constructs a view of the array's memory with a different data-type. This can cause a reinterpretation of the bytes of memory.

``a.view(ndarray_subclass)`` or ``a.view(type=ndarray_subclass)`` just returns an instance of ``ndarray_subclass`` that looks at the same array (same shape, dtype, etc.) This does not cause a reinterpretation of the memory.

For ``a.view(some_dtype)`` , if ``some_dtype`` has a different number of bytes per entry than the previous dtype (for example, converting a regular array to a structured array), then the behavior of the view cannot be predicted just from the superficial appearance of ``a`` (shown by ``print(a)``). It also depends on exactly how ``a`` is stored in memory. Therefore if ``a`` is C-ordered versus fortran-ordered, versus defined as a slice or transpose, etc., the view may give different results.

Examples

>>> x = np.array([(1, 2)], dtype=[('a', np.int8), ('b', np.int8)])

Viewing array data using a different type and dtype:

```
>>> y = x.view(dtype=np.int16, type=np.matrix)
>>> y
matrix([[513]], dtype=int16)
>>> print type(y)
<class 'numpy.matrixlib.defmatrix.matrix'>
```

Creating a view on a structured array so it can be used in calculations

```
>>> x = np.array([(1, 2),(3,4)], dtype=[('a', np.int8), ('b', np.int8)])
>>> xv = x.view(dtype=np.int8).reshape(-1,2)
>>> xv
array([[1, 2],
       [3, 4]], dtype=int8)
>>> xv.mean(0)
array([ 2.,  3.])
```

Making changes to the view changes the underlying array

```
>>> xv[0,1] = 20
>>> print x
[(1, 20) (3, 4)]
```

Using a view to convert an array to a record array:

```
>>> z = x.view(np.recarray)
>>> z.a
array([1], dtype=int8)
```

Views share data:

```
>>> x[0] = (9, 10)
>>> z[0]
(9, 10)
```

Views that change the dtype size (bytes per entry) should normally be avoided on arrays defined by slices, transposes, fortran-ordering, etc.:

```
>>> x = np.array([[1,2,3],[4,5,6]], dtype=np.int16)
>>> y = x[:, 0:2]
>>> y
```



```

array([[1, 2],
       [4, 5]], dtype=int16)
>>> y.view(dtype=[('width', np.int16), ('length', np.int16)])
Traceback (most recent call last):
  File "<stdin>", line 1, in <module>
ValueError: new type not compatible with array.
>>> z = y.copy()
>>> z.view(dtype=[('width', np.int16), ('length', np.int16)])
array([[1, 2],
       [4, 5]], dtype=[('width', '<i2'), ('length', '<i2')])

```

Data descriptors inherited from [numpy.ndarray](#):

T

Same as `transpose()`, except that `self` is returned if `self.ndim < 2`.

Examples

```

-----
>>> x = np.array([[1.,2.],[3.,4.]])
>>> x
array([[ 1.,  2.],
       [ 3.,  4.]])
>>> x.T
array([[ 1.,  3.],
       [ 2.,  4.]])
>>> x = np.array([1.,2.,3.,4.])
>>> x
array([ 1.,  2.,  3.,  4.])
>>> x.T
array([ 1.,  2.,  3.,  4.])

```

__array_interface__

Array protocol: Python side.

__array_priority__

Array priority.

__array_struct__

Array protocol: C-struct side.

base

Base object if memory is from some other object.

Examples

The base of an array that owns its memory is `None`:

```

>>> x = np.array([1,2,3,4])
>>> x.base is None
True

```

Slicing creates a view, whose memory is shared with `x`:

```

>>> y = x[2:]
>>> y.base is x
True

```

ctypes

An object to simplify the interaction of the array with the `ctypes` module.

This attribute creates an object that makes it easier to use arrays when calling shared libraries with the `ctypes` module. The returned object has, among others, `data`, `shape`, and `strides` attributes (see Notes below) which themselves return `ctypes` objects that can be used as arguments to a shared library.

Parameters

None

Returns

c : Python object

Possessing attributes data, shape, strides, etc.

See Also

numpy.ctypeslib

Notes

Below are the public attributes of this object which were documented in "Guide to NumPy" (we have omitted undocumented public attributes, as well as documented private attributes):

- * data: A pointer to the memory area of the array as a Python integer. This memory area may contain data that is not aligned, or not in correct byte-order. The memory area may not even be writeable. The array flags and data-type of this array should be respected when passing this attribute to arbitrary C-code to avoid trouble that can include Python crashing. User Beware! The value of this attribute is exactly the same as `self._array_interface_['data'][0]`.
- * shape (`c_intp*self.ndim`): A ctypes array of length `self.ndim` where the basetype is the C-integer corresponding to `dtype('p')` on this platform. This base-type could be `c_int`, `c_long`, or `c_longlong` depending on the platform. The `c_intp` type is defined accordingly in `numpy.ctypeslib`. The ctypes array contains the shape of the underlying array.
- * strides (`c_intp*self.ndim`): A ctypes array of length `self.ndim` where the basetype is the same as for the shape attribute. This ctypes array contains the strides information from the underlying array. This strides information is important for showing how many bytes must be jumped to get to the next element in the array.
- * data_as(obj): Return the data pointer cast to a particular c-types object. For example, calling `self._as_parameter_` is equivalent to `data_as(ctypes.c_void_p)`. Perhaps you want to use the data as a pointer to a ctypes array of floating-point data: `data_as(ctypes.POINTER(ctypes.c_double))`.
- * shape_as(obj): Return the shape tuple as an array of some other c-types type. For example: `shape_as(ctypes.c_short)`.
- * strides_as(obj): Return the strides tuple as an array of some other c-types type. For example: `strides_as(ctypes.c_longlong)`.

Be careful using the ctypes attribute - especially on temporary arrays or arrays constructed on the fly. For example, calling `((a+b).ctypes.data_as(ctypes.c_void_p))` returns a pointer to memory that is invalid because the array created as `(a+b)` is deallocated before the next Python statement. You can avoid this problem using either `c=a+b` or `ct=(a+b).ctypes`. In the latter case, `ct` will hold a reference to the array until `ct` is deleted or re-assigned.

If the ctypes module is not available, then the ctypes attribute of array objects still returns something useful, but ctypes objects are not returned and errors may be raised instead. In particular, the object will still have the `as_parameter` attribute which will return an integer equal to the data attribute.

Examples

```
>>> import ctypes
>>> x
array([[0, 1],
       [2, 3]])
>>> x.ctypes.data
30439712
>>> x.ctypes.data_as(ctypes.POINTER(ctypes.c_long))
<ctypes.LP_c_long object at 0x01F01300>
>>> x.ctypes.data_as(ctypes.POINTER(ctypes.c_long)).contents
c_long(0)
```

```
>>> x.ctypes.data_as(ctypes.POINTER(ctypes.c_longlong)).contents
c_longlong(4294967296L)
>>> x.ctypes.shape
<numpy.core._internal.c_long_Array_2 object at 0x01FFD580>
>>> x.ctypes.shape_as(ctypes.c_long)
<numpy.core._internal.c_long_Array_2 object at 0x01FCE620>
>>> x.ctypes.strides
<numpy.core._internal.c_long_Array_2 object at 0x01FCE620>
>>> x.ctypes.strides_as(ctypes.c_longlong)
<numpy.core._internal.c_longlong_Array_2 object at 0x01F01300>
```

data

Python buffer object pointing to the start of the array's data.

dtype

Data-type of the array's elements.

Parameters

None

Returns

d : numpy dtype object

See Also

numpy.dtype

Examples

```
>>> x
array([[0, 1],
       [2, 3]])
>>> x.dtype
dtype('int32')
>>> type(x.dtype)
<type 'numpy.dtype'>
```

flags

Information about the memory layout of the array.

Attributes

C_CONTIGUOUS (C)

The data is in a single, C-style contiguous segment.

F_CONTIGUOUS (F)

The data is in a single, Fortran-style contiguous segment.

OWNDATA (O)

The array owns the memory it uses or borrows it from another object.

WRITEABLE (W)

The data area can be written to. Setting this to False locks the data, making it read-only. A view (slice, etc.) inherits WRITEABLE from its base array at creation time, but a view of a writeable array may be subsequently locked while the base array remains writeable. (The opposite is not true, in that a view of a locked array may not be made writeable. However, currently, locking a base object does not lock any views that already reference it, so under that circumstance it is possible to alter the contents of a locked array via a previously created writeable view onto it.) Attempting to change a non-writeable array raises a RuntimeError exception.

ALIGNED (A)

The data and all elements are aligned appropriately for the hardware.

UPDATEIFCOPY (U)

This array is a copy of some other array. When this array is deallocated, the base array will be updated with the contents of this array.

FNC

F_CONTIGUOUS and not C_CONTIGUOUS.

FORC

F_CONTIGUOUS or C_CONTIGUOUS (one-segment test).

BEHAVED (B)

ALIGNED and WRITEABLE.

CARRAY (CA)
BEHAVED and C_CONTIGUOUS.
FARRAY (FA)
BEHAVED and F_CONTIGUOUS and not C_CONTIGUOUS.

Notes

The ``flags`` object can be accessed dictionary-like (as in ``a.flags['WRITEABLE']``), or by using lowercased attribute names (as in ``a.flags.writeable``). Short flag names are only supported in dictionary access.

Only the `UPDATEIFCOPY`, `WRITEABLE`, and `ALIGNED` flags can be changed by the user, via direct assignment to the attribute or dictionary entry, or by calling ``ndarray.setflags``.

The array flags cannot be set arbitrarily:

- `UPDATEIFCOPY` can only be set ``False``.
- `ALIGNED` can only be set ``True`` if the data is truly aligned.
- `WRITEABLE` can only be set ``True`` if the array owns its own memory or the ultimate owner of the memory exposes a writeable buffer interface or is a string.

Arrays can be both C-style and Fortran-style contiguous simultaneously. This is clear for 1-dimensional arrays, but can also be true for higher dimensional arrays.

Even for contiguous arrays a stride for a given dimension ``arr.strides[dim]`` may be *arbitrary* if ``arr.shape[dim] == 1`` or the array has no elements. It does *not* generally hold that ``self.strides[-1] == self.itemsize`` for C-style contiguous arrays or ``self.strides[0] == self.itemsize`` for Fortran-style contiguous arrays is true.

flat

A 1-D iterator over the array.

This is a ``numpy.flatiter`` instance, which acts similarly to, but is not a subclass of, Python's built-in iterator object.

See Also

`flatten` : Return a copy of the array collapsed into one dimension.

`flatiter`

Examples

```
>>> x = np.arange(1, 7).reshape(2, 3)
>>> x
array([[1, 2, 3],
       [4, 5, 6]])
>>> x.flat[3]
4
>>> x.T
array([[1, 4],
       [2, 5],
       [3, 6]])
>>> x.T.flat[3]
5
>>> type(x.flat)
<type 'numpy.flatiter'>
```

An assignment example:

```
>>> x.flat = 3; x
array([[3, 3, 3],
       [3, 3, 3]])
>>> x.flat[[1,4]] = 1; x
array([[3, 1, 3],
       [3, 1, 3]])
```

imag

The imaginary part of the array.

Examples

```
>>> x = np.sqrt([1+0j, 0+1j])
>>> x.imag
array([ 0.          ,  0.70710678])
>>> x.imag.dtype
dtype('float64')
```

itemsize

Length of one array element in bytes.

Examples

```
>>> x = np.array([1,2,3], dtype=np.float64)
>>> x.itemsize
8
>>> x = np.array([1,2,3], dtype=np.complex128)
>>> x.itemsize
16
```

nbytes

Total bytes consumed by the elements of the array.

Notes

Does not include memory consumed by non-element attributes of the array object.

Examples

```
>>> x = np.zeros((3,5,2), dtype=np.complex128)
>>> x.nbytes
480
>>> np.prod(x.shape) * x.itemsize
480
```

ndim

Number of array dimensions.

Examples

```
>>> x = np.array([1, 2, 3])
>>> x.ndim
1
>>> y = np.zeros((2, 3, 4))
>>> y.ndim
3
```

real

The real part of the array.

Examples

```
>>> x = np.sqrt([1+0j, 0+1j])
>>> x.real
array([ 1.          ,  0.70710678])
>>> x.real.dtype
dtype('float64')
```

See Also

numpy.real : equivalent function

shape

Tuple of array dimensions.

Notes

May be used to "reshape" the array, as long as this would not require a change in the total number of elements

Examples

```

-----
>>> x = np.array([1, 2, 3, 4])
>>> x.shape
(4,)
>>> y = np.zeros((2, 3, 4))
>>> y.shape
(2, 3, 4)
>>> y.shape = (3, 8)
>>> y
array([[ 0.,  0.,  0.,  0.,  0.,  0.,  0.,  0.],
       [ 0.,  0.,  0.,  0.,  0.,  0.,  0.,  0.],
       [ 0.,  0.,  0.,  0.,  0.,  0.,  0.,  0.]])
>>> y.shape = (3, 6)
Traceback (most recent call last):
  File "<stdin>", line 1, in <module>
ValueError: total size of new array must be unchanged

```

size

Number of elements in the array.

Equivalent to `np.prod(a.shape)`, i.e., the product of the array's dimensions.

Examples

```

-----
>>> x = np.zeros((3, 5, 2), dtype=np.complex128)
>>> x.size
30
>>> np.prod(x.shape)
30

```

strides

Tuple of bytes to step in each dimension when traversing an array.

The byte offset of element `(i[0], i[1], ..., i[n])` in an array `a` is::

```
offset = sum(np.array(i) * a.strides)
```

A more detailed explanation of strides can be found in the "ndarray.rst" file in the NumPy reference guide.

Notes

Imagine an array of 32-bit integers (each 4 bytes)::

```
x = np.array([[0, 1, 2, 3, 4],
              [5, 6, 7, 8, 9]], dtype=np.int32)
```

This array is stored in memory as 40 bytes, one after the other (known as a contiguous block of memory). The strides of an array tell us how many bytes we have to skip in memory to move to the next position along a certain axis. For example, we have to skip 4 bytes (1 value) to move to the next column, but 20 bytes (5 values) to get to the same position in the next row. As such, the strides for the array `x` will be `(20, 4)`.

See Also

`numpy.lib.stride_tricks.as_strided`

Examples

```

>>> y = np.reshape(np.arange(2*3*4), (2,3,4))
>>> y
array([[[ 0,  1,  2,  3],
        [ 4,  5,  6,  7],
        [ 8,  9, 10, 11]],
       [[12, 13, 14, 15],
        [16, 17, 18, 19],
        [20, 21, 22, 23]]])
>>> y.strides
(48, 16, 4)
>>> y[1,1,1]

```

```
17
>>> offset=sum(y.strides * np.array((1,1,1)))
>>> offset/y.itemsize
17

>>> x = np.reshape(np.arange(5*6*7*8), (5,6,7,8)).transpose(2,3,1,0)
>>> x.strides
(32, 4, 224, 1344)
>>> i = np.array([3,5,2,2])
>>> offset = sum(i * x.strides)
>>> x[3,5,2,2]
813
>>> offset / x.itemsize
813
```

Functions

loadbpf(fname)

Loads biosignal data.

Parameters

fname : string
input filename

Returns

bpd : [bpcarray](#)
output biosignal array

See Also

[np.loadtxt](#)

Notes

Example

References

.. [1]

<http://>

```
.. module:: resp
   :platform: Unix, Windows
   :synopsis: This module provides various functions to handle RESP signals.

.. moduleauthor:: Filipe Canento
```

Modules

[numpy](#)
[scipy](#)
[unittest](#)
[pylab](#)
[sys](#)

Functions

features(Signal=None, SamplingRate=1000.0, Filter={})
Retrieves relevant RESP signal features.

Kwargs:

- Signal (array): input signal.
- SamplingRate (float): Sampling frequency (Hz).
- Filter (dict): filter parameters.

Kwrvals:

- Signal (array): output filtered signal.
- FR (array): instant respiratory frequency (Hz)
- ZC (array): zero crossings indexes
- mean (float): mean
- std (float): standard deviation
- var (float): variance
- skew (ndarray): skewness
- kurtosis (array): kurtosis
- ad (float): absolute deviation

Configurable fields:
{"name": "resp.features", "config": {"SamplingRate": "1000."}, "inputs": ["Signal", "Filter"], "outputs": ["Signal", "FR", "ZC", "mean", "std", "var", "skew", "kurtosis", "ad"]}

See Also:

- filt
- tls.statsf

Notes:

Example:

References:

- .. [1]

filt(Signal=None, SamplingRate=1000.0, UpperCutoff=0.35, LowerCutoff=0.1, Order=2.0)
Filters an input RESP signal.

If only input signal is provide, it returns the filtered RESP signal assuming a 1000Hz sampling frequency and ...

Kwargs:

- Signal (array): input signal.
- SamplingRate (float): Sampling frequency (Hz).
- UpperCutoff (float): Low-pass filter cutoff frequency (Hz).
- LowerCutoff (float): High-pass filter cutoff frequency (Hz).
- Order (int): Filter order.

Kwrvals:

- Signal (array): output filtered signal.

Configurable fields:
{"name": "resp.filt", "config": {"UpperCutoff": "0.35", "SamplingRate": "1000.", "LowerCutoff": "0.1", "Order": "2."}, "inputs": ["Signal"], "outputs": ["Signal"]}

See Also:

- flt.zpdf

Notes:

Example:
Signal = load(...)
SamplingRate = ...
res = [filt](#)(Signal=Signal, SamplingRate=SamplingRate)
plot(res['Signal'])

References:
.. [1]

resp(Signal=None, SamplingRate=1000.0, Filter=())
Respiratory signal information.

Kwargs:
Signal (array): input signal.

SamplingRate (float): Sampling frequency (Hz).

Filter (dict): filter parameters.

Kwrvals:
Signal (array): output filtered signal.

FR (array): instant respiratory frequency (Hz)

ZC (array): zero crossings indexes

Configurable fields:{"name": "resp.resp", "config": {"SamplingRate": "1000.0"}, "inputs": ["Signal", "Filter"], "outputs": ["Signal", "FR", ']

See Also:
[filt](#)

Notes:

Example:

References:
.. [1]

This module provides various functions to ...

Functions:

[step\(\)](#)
[threshold\(\)](#)
[match\(\)](#)
[pair\(\)](#)

Modules

[pylab](#)

Functions

match(Signal=None, Window=None)

.

Parameters

Signal : ndarray
Input signal data.
Window :

Returns

kwrvals : dict
A keyworded return values dict is returned with the following keys:
Event :

See Also

pl.convolve

Notes

Example

pair(x, y)

.

Parameters

x : array

y : array

Returns

x : array

y : array

See Also

Notes

Example

step(Signal=None, Shift='>')

Find the indexes where a dirac, unit step or unit pulse function in the input signal rises or falls.

This implementation is based on the discrete difference of the input signal.

Parameters

Signal : ndarray

Input signal data with a dirac, unit step or unit pulse function.

Shift : str

Direction of detection, '>' for rise or '<' for fall.

Default: '>'

Returns

kwrvls : dict

A keyworded return values dict is returned with the following keys:

Event : ndarray

The indexes within the input signal 'Signal' where the 'Shift' was detected.

See Also

sync.threshold

sync.match

Example

```
x = zeros(6)
x[2:4] = 1
plot(x)
vlines(step(x)['Events'],min(x),max(x),'r','dashed')
vlines(step(x,'<')['Events'],min(x),max(x),'g','dashed')
legend(('unit pulse','rise','fall'))
```

References

.. [1] Wikipedia, "Dirac Delta Function".

http://en.wikipedia.org/wiki/Dirac_delta_function

.. [2] Wikipedia, "Heaviside Step Function".

http://en.wikipedia.org/wiki/Heaviside_step_function

.. [3] Wikipedia, "Unit Pulse Function".

http://en.wikipedia.org/wiki/Rectangular_function

threshold(Signal=None, Threshold=0.1, Shift='>')

Find the indexes where a the input signal rises above or falls bellow a given threshold.

Parameters

Signal : ndarray

Input signal data.

Threshold : int, float

Detection threshold. If this is an 'int' greater or equal than one it will be directly used as the threshold. If this is a 'float', the threshold is computed as the 'Threshold' percentage of the signal span.

Default: .1
Shift : str
Direction of detection, `>` for rise or `<` for fall.
Default: `>`

Returns

kwargs : dict
A keyworded return values dict is returned with the following keys:
Event : ndarray
The indexes within the input signal `Signal` where a `Shift` with respect to `Threshold` was detected.

See Also

sync.step
sync.match

Notes

This method is primarily designed to detect shifts in the input signal with respect to a threshold computed as a percentage of the signal span.

Example

```
x = zeros(6)
x[2:4] = 1
plot(x)
th=.5
plot(ones(len(x))*th,'k--')
vlines(threshold(x,th)['Event'],min(x),max(x),'r','dashed')
vlines(threshold(x,th,'<')['Event'],min(x),max(x),'g','dashed')
legend(('unit pulse','rise', 'fall'))
```

This module provides various auxiliary functions.

Functions:

[statsf\(\)](#)
[zerocross\(\)](#)
[scale\(\)](#)
[walktree\(\)](#)
[fileparts\(\)](#)
[fullfile\(\)](#)
[format_num\(\)](#)
[get_max_width\(\)](#)
[pprint_table\(\)](#)

Modules

[locale](#)
[numpy](#)

[os](#)
[pylab](#)

[sys](#)

Functions

biplot(m)

fileparts(file)

#-----

format_num(num)

Format a number according to given places.
Adds commas, etc. Will truncate floats into ints!

fullfile(*args)

#-----

get_max_width(table, index)

Get the maximum width of the given column index

pprint_table(out, table)

Prints out a table of data, padded for alignment
@param out: Output stream (file-like object)
@param table: The table to print. A list of lists.
Each row must have the same number of columns.

scale(x)

#-----

statsf(Signal=None)

Determine various statistical features of a given input signal.

Parameters

Signal : array
input signal

Returns

kwrvls : dict
A keyworded return values dict is returned with the following keys:

- 'mean' : float
mean
- 'std' : float
standard deviation
- 'var' : float
variance
- 'skew' : ndarray
skewness
- 'kurtosis' : array
kurtosis
- 'ad' : float
absolute deviation

See Also

- np.mean
- np.std
- np.var
- np.skew
- scipy.stats.skew
- scipy.stats.kurtosis

Notes

Example

References

- .. [1]
- .. [2]
- <http://>

walktree(top, callback=<function disp>)

recursively descend the directory tree rooted at top,
calling the callback function for each regular file

EXTRACTED FROM:

Python Library

<http://docs.python.org/library/stat.html>

zerocross(Signal=None)

Determine the indexes where the signal crosses zero.

Parameters

signal : array
input signal

Returns

kwrvals : dict

A keyworded return values dict is returned with the following keys:
 'zc' : array
 indexes of the zero crossings.

See Also

Notes

Example

References

.. [1]
.. [2]
<http://>

Data

SF_APPEND = 262144
SF_ARCHIVED = 65536
SF_IMMUTABLE = 131072
SF_NOUNLINK = 1048576
SF_SNAPSHOT = 2097152
ST_ETIME = 7
ST_CTIME = 9
ST_DEV = 2
ST_GID = 5
ST_INO = 1
ST_MODE = 0
ST_MTIME = 8
ST_NLINK = 3
ST_SIZE = 6
ST_UID = 4
S_ENFMT = 1024
S_IXEXEC = 64
S_IFBLK = 24576
S_IFCHR = 8192
S_IFDIR = 16384
S_IFIFO = 4096
S_IFLNK = 40960
S_IFREG = 32768
S_IFSOCK = 49152

S_IREAD = 256
S_IRGRP = 32
S_IROTH = 4
S_IRUSR = 256
S_IRWXG = 56
S_IRWXO = 7
S_IRWXU = 448
S_ISGID = 1024
S_ISUID = 2048
S_ISVTX = 512
S_IWGRP = 16
S_IWOTH = 2
S_IWRITE = 128
S_IWUSR = 128
S_IXGRP = 8
S_IXOTH = 1
S_IXUSR = 64
UF_APPEND = 4
UF_COMPRESSED = 32
UF_HIDDEN = 32768
UF_IMMUTABLE = 2
UF_NODUMP = 1
UF_NOUNLINK = 16
UF_OPAQUE = 8

bspyclone.bvp

[index](#)

[d:\work\productioncode\clones\biosppy\bspyclone\bvp__init__.py](#)

```
# import bvp
# import models
```

Package Contents

[bvp](#)

[models](#)

```
.. module:: h5db
   :platform: Unix, Windows
   :synopsis: This module provides a wrapper to the HDF5 file format, adapting it to store biosignals according to the BioMESH specification at http://cam

.. moduleauthor:: Carlos Carreiras
```

Modules

[h5py](#) [json](#)

Classes

[IOController](#)
[hdf](#)
[meta](#)
[realTime](#)

class IOController

Methods defined here:

```
__init__(self, fileConfig)

put(self, data)

start(self)

stop(self)
```

class hdf

Wrapper class to operate on HDF5 records according to the BioMESH specification.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
 .. [1]

Methods defined here:

```
__enter__(self)
    __enter__ Method for 'with' statement.

    Kwargs:
        None

    Kwrvals:
        None

    See Also:

    Notes:

    Example:

    References:
        .. [1]

__exit__(self, exc_type, exc_value, traceback)
    __exit__ Method for 'with' statement.

    Kwargs:
        None

    Kwrvals:
        None
```

See Also:

Notes:

Example:

References:
.. [1]

__init__(self, filePath=None, mode='a')
Open the HDF5 record.

Kwargs:
filePath (str): Path to HDF5 file.

mode (str): File access mode. Available modes:
 'r+': Read/write, file must exist
 'r': Read only, file must exist
 'w': Create file, truncate if exists
 'w-': Create file, fail if exists
 'a': Read/write if exists, create otherwise
 Default: 'a'.

Kwrvals:

See Also:

Notes:

Example:
fid = [hdf](#)('record.hdf5', 'a')

References:
.. [1]

addEvent(self, eventType="", timeStamps=None, values=None, mdata=None, eventName=None, compress=False)
Method to add asynchronous data (events) to the HDF5 record.

Kwargs:
eventType (str): Type of the events to add. Default: ''.

timeStamps (array): Array of time stamps. Default: [].

values (array): Array with data for each time stamp. Default: [].

mdata (dict): Dictionary object with metadata about the events. Default: {}.

eventName (str): Name of the group to be created.

compress (bool): Flag to compress the data (GZIP). Default: False.

Kwrvals:

See Also:

Notes:

Example:
fid.[addEvent](#)('/test', [0, 1, 2], [[1, 2], [3, 4], [5, 6]], {'comments': 'test event'}, 'event0')

References:
.. [1]

addInfo(self, header={})
Method to add or overwrite the basic information (header) of the HDF5 record.

Kwargs:
header (dict): Dictionary (JSON) object with the information. Default: {}.

Kwrvals:

See Also:

Notes:

Example:
fid.[addInfo](#)({'name': 'record'})

References:
.. [1]

addSignal(self, signalType="", signal=None, mdata=None, dataName=None, compress=False)
Method to add a signal (synchronous data) to the HDF5 record.

Kwargs:

- signalType (str): Type of the signal to add. Default: ''.
- signal (array): Array with the signal to add.
- mdata (dict): Dictionary object with metadata about the data. Default: {}.
- dataName (str): Name of the dataset to be created.
- compress (bool): Flag to compress the data (GZIP). Default: False.

Kwrvals:

See Also:

Notes:

Example:

```
fid.addSignal('/test', [0, 1, 2], {'comments': '/test signal'}, 'signal0')
```

References:

.. [1]

addSignalRT(self, signalType="", mdata=None, dataName=None, blockShape=None, axis=0, dtype='f8', compress=False)
Method to add a signal (synchronous data) to the HDF5 record in real time.

Kwargs:

- signalType (str): Type of the signal to add. Default: ''.
- mdata (dict): Dictionary object with metadata about the data. Default: {}.
- dataName (str): Name of the dataset to be created.
- blockShape (tuple): Shape of the signal blocks.
- axis (int): Direction on which the expansion of the dataset is made. Default: 0.
- dtype (str): Data type of the signal to add (supports numpy data types).
- compress (bool): Flag to compress the data (GZIP). Default: False.

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

close(self)
Method to close the HDF5 record.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

```
fid.close()
```

References:

.. [1]

delEvent(self, eventType="", eventName=None)
Method to delete asynchronous data (events) from the HDF5 record. The record is marked for repackaging.

Kwargs:

- eventType (str): Type of the desired event. Default: ''.
- eventName (str): Name of the dataset to retrieve.

Kwrvals:

See Also:

Notes:

Example:
fid.[delEvent](#)('/test', 'event0')

References:
.. [1]

delEventType(self, eventType="")

Method to delete a type of asynchronous data from the HDF5 record. The record is marked for repackaging.

Kwargs:
 eventType (str): Type of the desired event. Default: ''.

Kwrvals:

See Also:

Notes:

Example:
fid.[delEventType](#)('/test')

References:
.. [1]

delSignal(self, signalType="", dataName=None)

Method to delete a signal (synchronous data) from the HDF5 record. The record is marked for repackaging.

Kwargs:
 signalType (str): Type of the desired signal. Default: ''.
 dataName (str): Name of the dataset to retrieve.

Kwrvals:

See Also:

Notes:

Example:
fid.[delSignal](#)('/test', 'data0')

References:
.. [1]

delSignalType(self, signalType="")

Method to delete a type of signals from the HDF5 record. The record is marked for repackaging.

Kwargs:
 signalType (str): Signal type to delete. Default: ''.

Kwrvals:

See Also:

Notes:

Example:
fid.[delSignalType](#)('/test')

References:
.. [1]

getEvent(self, eventType="", eventName=None)

Method to retrieve asynchronous data(events) from the HDF5 record.

Kwargs:
 eventType (str): Type of the desired event. Default: ''.
 eventName (str): Name of the dataset to retrieve.

Kwrvals:
 timeStamps (array): Array of time stamps.
 values (array): Array with data for each time stamp.
 mdata (dict): Dictionary object with metadata about the events.

See Also:

Notes:

Example:
out = fid.[getEvent](#)('/test', 'event0')
timeStamps = out['timeStamps']

```
values = out['values']
metadata = out['mdata']
```

References:
.. [1]

getEventInfo(self, eventWeg=None)

Method to retrieve the metadata of asynchronous.

Kwargs:
eventWeg (str): Path to the group.

Kwrvals:
mdata (dict): Dictionary with the desired information.

See Also:

Notes:

Example:
out = fid.[getEventInfo](#)('/test/event0')
metadata = out['mdata']

References:
.. [1]

getInfo(self)

Method to retrieve the basic information (header) of the HDF5 record.

Kwargs:

Kwrvals:
header (dict): Dictionary object with the header information.

See Also:

Notes:

Example:
header = fid.get()['header']

References:
.. [1]

getRepack(self)

Get the repack flag.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:
repack = fid.[getRepack](#)()

References:
.. [1]

getSignal(self, signalType="", dataName=None)

Method to retrieve a signal (synchronous data) from the HDF5 record.

Kwargs:
signalType (str): Type of the desired signal. Default: ''.
dataName (str): Name of the dataset to retrieve.

Kwrvals:
signal (array): Array with the signals.
mdata (dict): Dictionary object with metadata about the signals.

See Also:

Notes:

Example:
out = fid.[getSignal](#)('/test', 'signal0')
signal = out['signal']
metadata = out['mdata']

References:
.. [1]

getSignalInfo(self, dataWeg=None)

Method to retrieve the metadata of a signal.

Kwargs:

dataWeg (str): Path to the dataset.

Kwrvals:

mdata (dict): Dictionary with the desired information.

See Also:

Notes:

Example:

```
out = fid.getSignalInfo('/test/data0')
metadata = out['mdata']
```

References:

.. [1]

getSignalSet(self, signalType="", dataName=None)

Method to retrieve a signal (synchronous data) from the HDF5 record.

Kwargs:

signalType (str): Type of the desired signal. Default: ''.

dataName (str): Name of the dataset to retrieve.

Kwrvals:

signal (array): Array with the signals.

mdata (dict): Dictionary object with metadata about the signals.

See Also:

Notes:

Example:

```
out = fid.getSignal('/test', 'signal0')
signal = out['signal']
metadata = out['mdata']
```

References:

.. [1]

listEvents(self, eventType="")

Method to list the events (asynchronous data) belonging to a given type.

Kwargs:

eventType (str): Type of the desired signal. Default: ''.

Kwrvals:

eventsList (lit): List of the events of the given type.

See Also:

Notes:

Example:

```
fid.listEvents('/test')
```

References:

.. [1]

listSignals(self, signalType="")

Method to list the signals (synchronous data) belonging to a given type.

Kwargs:

signalType (str): Type of the desired signal. Default: ''.

Kwrvals:

signalsList (lit): List of the signals of the given type.

See Also:

Notes:

Example:

```
fid.listSignals('/test')
```

References:

.. [1]

setRepack(self)

Set flag to repack.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:
fid.[setRepack\(\)](#)

References:
.. [1]

unsetRepack(self)
Unset flag to repack.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:
fid.[unsetRepack\(\)](#)

References:
.. [1]

class meta

Wrapper class to store experiments and subjects on HDF5.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

Methods defined here:

__enter__(self)
__enter__ Method for 'with' statement.

Kwargs:
None

Kwrvals:
None

See Also:

Notes:

Example:

References:
.. [1]

__exit__(self, exc_type, exc_value, traceback)
__exit__ Method for 'with' statement.

Kwargs:
None

Kwrvals:
None

See Also:

Notes:

Example:

References:

.. [1]

__init__(self, filePath=None, mode='a')
Open the HDF5 file.

Kwargs:

filePath (str): Path to HDF5 file.

mode (str): File access mode. Available modes:

'r+': Read/write, file must exist

'r': Read only, file must exist

'w': Create file, truncate if exists

'w-': Create file, fail if exists

'a': Read/write if exists, create otherwise

Default: 'a'.

Kwrvals:

See Also:

Notes:

Example:

fid = [meta](#)('expsub.hdf5', 'a')

References:

.. [1]

addExperiment(self, experiment={})
Method to add an experiment to the file.

Kwargs:

experiment (dict): Dictionary with the experiment information. Default: {}.

Kwrvals:

See Also:

Notes:

The experiment must have a 'name' key.

Example:

fid.[addExperiment](#)({ 'name': 'experiment', 'comments': 'Hello world.'})

References:

.. [1]

addSubject(self, subject={})
Method to add a subject to the file.

Kwargs:

subject (dict): Dictionary with the subject information. Default: {}.

Kwrvals:

See Also:

Notes:

The subject must have a '_id' key.

Example:

fid.[addSubject](#)({ '_id': 0, 'name': 'subject'})

References:

.. [1]

close(self)
Method to close the HDF5 file.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

`fid.close()`

References:
.. [1]

getExperiment(self, experimentName=None)

Method to get the information about an experiment.

Kwargs:
 experimentName (str): The name of the experiment.

Kwrvals:
 experiment (dict): Dictionary with the experiment information.

See Also:

Notes:

Example:
 experiment = fid.[getExperiment](#)('experiment')['experiment']

References:
.. [1]

getSubject(self, subjectId=None)

Method to get the information about a subject.

Kwargs:
 subjectId (int): The ID of the subject.

Kwrvals:
 subject (dict): Dictionary with the subject information.

See Also:

Notes:

Example:
 subject = fid.[getSubject](#)(0)['subject']

References:
.. [1]

listExperiments(self)

Method to list all the experiments in the file.

Kwargs:

Kwrvals:
 expList (list): List with the experiments.

See Also:

Notes:

Example:
 expList = fid.[listExperiments](#)()['expList']

References:
.. [1]

listSubjects(self)

Method to list all the subjects in the file.

Kwargs:

Kwrvals:
 subList (list): List with the subjects.

See Also:

Notes:

Example:
 subList = fid.[listSubjects](#)()['subList']

References:
.. [1]

setDB(self, dbName=None)

Method to set the DB the file belongs to.

Kwargs:
 dbName (str): Name of the database.

Kwrvals:

See Also:

Notes:

Example:

```
fid.setDB('database')
```

References:

.. [1]

updateExperiment(self, experimentName=None, info={})

Method to update an experiment's information.

Kwargs:

experimentName (str): Name of the experiment.

info (dict): Dictionary with the information to update.

Kwrvals:

See Also:

Notes:

Example:

```
fid.updateExperiment('experiment', {'new': 'field'})
```

References:

.. [1]

updateSubject(self, subjectId=None, info={})

Method to update a subject's information.

Kwargs:

subjectId (int): The ID of the subject.

info (dict): Dictionary with the information to update.

Kwrvals:

See Also:

Notes:

Example:

```
fid.updateSubject(0, {'new': 'field'})
```

References:

.. [1]

class realTime

Class to add signals in real time to an HDF5 file.

Kwargs:

Kwrvals:

See Also:

Notes:

Example:

References:

.. [1]

Methods defined here:

__init__(self, group, mdata, dataName, blockShape, axis, dtype, compress=False)
Initialize the new dataset.

Kwargs:

group (h5py.Group, h5py.File): File or group instance.

mdata (dict): Dictionary object with metadata about the data.

dataName (str): Name of the dataset to be created.

`blockShape (tuple)`: Shape of the signal blocks.

`axis (int)`: Direction on which the expansion of the dataset is made.

`dtype (str)`: Data type of the signal to add (supports numpy data types).

`compress (bool)`: Flag to compress the data (GZIP). Default: False.

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

put(self, data)
Add data.

Kwargs:
`data (array)`: Data to add.

Kwrvals:

See Also:

Notes:

Example:

References:
.. [1]

Functions

IOProcess(fileConfig, acquisition, queue)

latin1ToAscii(text)

This replaces UNICODE Latin-1 characters with something equivalent in 7-bit ASCII.

Kwargs:
`text (str)`: Text to convert.

Kwrvals:

See Also:

Notes:

Example:

```
text = 'João'
print text
ntext = latin1ToAscii(text)
print ntext
```

References:
.. [1] <http://stackoverflow.com/questions/930303/python-string-cleanup-manipulation-accented-characters>

bspyclone.database [index](#)
[d:\work\productioncode\clones\biosppy\bspyclone\database__init__.py](#)

```
# import biomes  
# import h5db
```

Package Contents

[HDF5DB](#)
[biomes](#)
[example](#)

[h5db](#)
[h5repack](#)
[mongoH5](#)

[old2biomes](#)
[philipsXML](#)
[physionet](#)

[repackScript](#)
[sandbox](#)
[syncdb](#)

bspyclone.ecg

[index](#)

[d:\work\productioncode\clones\biosppy\bspyclone\ecg__init__.py](#)

```
# import ecg
# import models
# import tools
```

Package Contents

[ecg](#)
[hrv](#)

[models](#)
[models_francis](#)

[models_francis_old](#)
[tools](#)

bspyclone.eda

[index](#)

[d:\work\productioncode\clones\biosppy\bspyclone\eda__init__.py](#)

```
# import eda
# import models
```

Package Contents

[eda](#)

[models](#)

bspyclone.eeg

[index](#)

[d:\work\productioncode\clones\biosppy\bspyclone\eeg__init__.py](#)

```
# import eeg
# reload(eeg)
# from eeg import *
```

Package Contents

[eeg](#)

bspyclone.emg

[index](#)

[d:\work\productioncode\clones\biosppy\bspyclone\emg__init__.py](#)

```
# import emg
```

Package Contents

[emg](#)

bspyclone

[index](#)

[d:\work\productioncode\clones\biosppy\bspyclone__init__.py](#)

BIOSPPY.

.. moduleauthor:: Filipe Canento, Hugo Silva, Andre Lourenco, Prof. Ana Fred

Package Contents

__apiBP	ecg (package)	peakd	sandbox
apiBP	eda (package)	plx	sync
bioplux	eeg (package)	plxv2	timing
biopluxwx	emg (package)	pywx	tools
bvp (package)	filt	resp (package)	
database (package)	hr	sand	

bspyclone.resp

[index](#)

[d:\work\productioncode\clones\biosppy\bspyclone\resp__init__.py](#)

```
# import resp
```

Package Contents

[resp](#)